

NūPCR™ Assay Protocol Guide

FOR RESEARCH USE ONLY

Introduction	3
Best Practices	6
Contents	7
Consumables and Equipment	8
Assay Workflow	9
Assay Protocol	10
Data Analysis	13
Technical Assistance	



This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE) OR ANY USE OF SUCH PRODUCT(S) OUTSIDE THE SCOPE OF THE EXPRESS WRITTEN LICENSES OR PERMISSIONS GRANTED BY ILLUMINA IN CONNECTION WITH CUSTOMER'S ACQUISITION OF SUCH PRODUCT(S).

FOR RESEARCH USE ONLY

© 2012 Illumina, Inc. All rights reserved.


Illumina, illuminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CPro, DASL, DesignStudio, Eco, GAIx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

Dye Compounds in this product are sold under license from Biosearch Technologies, Inc. and protected by U.S. and world-wide patents either issued or in application.

Introduction

This document describes the NuPCR™ Assay protocol used to amplify and detect any DNA target (genomic DNA (gDNA), complementary DNA (cDNA) or plasmid DNA). NuPCR Gene Expression assays are designed using Illumina® DesignStudio™ software and contain all the oligonucleotides required for performing nucleic acid enzyme, NuZyme™, based quantitative PCR (qPCR). The NuPCR Assay can be run on any real-time PCR instrument.

Figure 1 NuPCR Assay Concept



NuZymes are nucleic acid enzymes that recognize and assemble on target DNA sequences forming a catalytic complex. The enzymatic activity of the NuZyme cleaves a fluorescently-labeled universal substrate, producing a signal that can be detected by a real-time PCR instrument.

NuZymes are composed of two oligonucleotides that are partial enzymes or PartZymes™. Each PartZyme contains part of the catalytic core sequence of the NuZyme, flanked by a substrate arm sequence and a sensor arm sequence. Each

sensor arm binds with high specificity to its target sequence, while the substrate arms bind the fluorescently-labeled, quenched universal substrate.

PartZymes contain no catalytic activity individually or in the absence of the target sequence. The target sequence serves to bring the two partzymes in close proximity, which facilitates the formation of the catalytically active NuZyme. The active NuZyme can now cleave the universal substrate oligo, however the target sequence remains intact. NuZymes are multiple-turnover enzymes, cleaving multiple substrate oligos in succession.

During a NuPCR reaction, target sequences are amplified using 5' and 3' primers and a DNA polymerase. The sensor arms of the partzymes bind to their corresponding sequences in the target amplicon, forming the catalytic NuZyme. The universal substrate binds the NuZyme substrate arms and is cleaved, releasing the fluorophore from the quencher, resulting in increased fluorescence. The fluorescent signal allows for detection and quantification of target sequences in real time.

NuPCR assays are supplied as a 20X solution containing all the oligonucleotides required for performing NuZyme based qPCR. The NuPCR master mix is provided as a 2X formulation containing all the components required for real-time PCR except primers, universal substrate, and DNA template. The ROX passive reference dye can be used for instruments that require it.



NOTE

NuPCR uses a three-step cycling protocol. Please read the *Assay Protocol* on page 10 carefully to properly setup your PCR instrument with the correct cycling conditions and fluorescence data acquisition step.



NOTE

- The NuPCR Assay can be run on any real-time PCR instrument. It has been shown to work on the following:
- Eco™ Real-Time PCR System
- Applied Biosystems 7900HT
- Applied Biosystems 7500
- Applied Biosystems ViiA 7
- Bio-Rad CFX 384



NOTE

You can access DesignStudio on the Illumina website at <http://www.illumina.com/applications/designstudio.ilmn>. A MyIllumina account is required.

Multiplex PCR with NuPCR

Two-color multiplex real-time PCR can be performed with NuPCR. DesignStudio can design optimal primer and NuZyme sets for single tube multiplex assays for up to two targets. All sets are analyzed for oligo interactions and homology to other targets for specific and efficient amplification.

What's New

The following changes have been made in this guide revision:

- ▶ Added a note in the Introduction to emphasize the importance of proper programming of the PCR instrument.
- ▶ Moved the workflow diagram to directly precede the protocol.
- ▶ Added a PCR instrument program diagram to the protocol.

Best Practices

When performing qPCR, you should always adhere to good molecular biology practices to prevent contamination, which can lead to false positive results.

- ▶ Wear a clean lab coat and clean gloves when preparing samples.
- ▶ Change gloves between handling different samples.
- ▶ Clean work surfaces thoroughly before and after the procedure.
- ▶ Maintain separate areas and segregate equipment and supplies for pre-PCR setup and post-PCR analysis.
 - Do not bring amplified PCR products into the pre-PCR area.
- ▶ Use aerosol-resistant filter tips.
- ▶ Pipette carefully to avoid spillage.
- ▶ Protect the NuPCR Gene Expression Assay mix from light. Excessive exposure to light may affect the fluorescent probe contained in the mix.
- ▶ Centrifuge the tube or plate before placing it in the real-time PCR instrument.

Contents

Check to ensure that you have at least one assay and one master mix identified in this section before proceeding.

Assays

Each NuPCR assay is shipped at room temperature. As soon as you receive it, store it at -15° to -25°C.

Table 1 NuPCR Assays

Description	Reactions	Catalog #
NuPCR Gene Expression Assays - FAM	250	EC-309-1001
NuPCR Gene Expression Assays - FAM	500	EC-309-1002
NuPCR Gene Expression Assays - FAM	2000	EC-309-1003
NuPCR Gene Expression Assays - HEX	250	EC-309-1011
NuPCR Gene Expression Assays - HEX	500	EC-309-1012
NuPCR Gene Expression Assays - HEX	2000	EC-309-1013
NuPCR Gene Expression Assays - ROX	250	EC-309-1021
NuPCR Gene Expression Assays - ROX	500	EC-309-1022
NuPCR Gene Expression Assays - ROX	2000	EC-309-1023
NuPCR Gene Expression Assays - Q670	250	EC-309-1031
NuPCR Gene Expression Assays - Q670	500	EC-309-1032
NuPCR Gene Expression Assays - Q670	2000	EC-309-1033

Master Mix

Each NuPCR master mix is shipped on dry ice. As soon as you receive it, store it at -15° to -25°C.

Table 2 NuPCR Master Mixes

Description	Reactions	Catalog #
NuPCR Gene Expression Master Mix	250	EC-305-1001
NuPCR Gene Expression Master Mix	500	EC-305-1002
NuPCR Gene Expression Master Mix	2000	EC-305-1003

Consumables and Equipment

Check to ensure that you have all of the necessary user-supplied consumables and equipment before proceeding to the NuPCR protocol.

Table 3 User-Supplied Consumables

Consumable	Supplier
2 µl barrier pipette tips	General lab supplier
2 µl multichannel pipettes	General lab supplier
2 µl single channel pipettes	General lab supplier
10 µl barrier pipette tips	General lab supplier
10 µl multichannel pipettes	General lab supplier
10 µl single channel pipettes	General lab supplier
Aerosol-resistant filter tips	General lab supplier
Optically clear adhesive seals	General lab supplier
PCR grade water	General lab supplier
PCR plates	General lab supplier
PCR tubes	General lab supplier
RNase/DNase zapper (to decontaminate surfaces)	General lab supplier


Table 4 User-Supplied Equipment

Equipment	Supplier
Real-time qPCR System	General lab supplier
Microplate centrifuge	General lab supplier
Vortexer	General lab supplier

Assay Workflow

The following figure illustrates the processes of the NuPCR Assay protocol to amplify and detect any DNA target.

Figure 2 NuPCR Assay Workflow



Assay Protocol

This protocol describes how to set up for a 20 μ l reaction volume. Reaction volumes can be scaled as needed for the real-time PCR instrument being used.

Illustrina-Supplied Consumables

- ▶ NuPCR Gene Expression Assays
- ▶ NuPCR Gene Expression Master Mix


User-Supplied Consumables

- ▶ DNA Template
- ▶ PCR Grade Water
- ▶ PCR Plate
- ▶ PCR Plate Adhesive Seal
- ▶ PCR Tube

Preparation

- ▶ Remove one each of the NuPCR Gene Expression Assays and NuPCR Gene Expression Master Mix tubes from -15° to -25°C storage and let stand to thaw at room temperature.
- ▶ Vortex the thawed NuPCR Gene Expression Assays and NuPCR Gene Expression Master Mix tubes to mix thoroughly.
- ▶ Pre-program the real-time PCR instrument thermal cycle settings as follows:
 - 95°C for 2 minutes
 - 40 cycles of:
 - 95°C for 15 seconds
 - 50°C for 30 seconds
 - 72°C for 30 seconds
- ▶ Program the real-time PCR instrument to acquire fluorescence data on each cycle of the 50°C , 30 second annealing step.

Figure 3 PCR Instrument Settings



Procedure

- Do one of the following:
 - For non-multiplex PCR**, prepare the following reaction mix in a new PCR tube. Vortex the tube to mix thoroughly. Multiply each volume by the number of samples being prepared. If you are preparing multiple samples, multiply each reagent volume by 1.2.

Reagent	Volume (μl)
NuPCR Gene Expression Master Mix	10
NuPCR Gene Expression Assays Mix	1
PCR Grade Water	enough to bring the total volume to 20 μl per sample
Total Volume per Sample	20 - DNA volume

- For multiplex PCR**, prepare the following reaction mix in a new PCR tube. Vortex the tube to mix thoroughly. Multiply each volume by the number of samples being prepared. If you are preparing multiple samples, multiply each reagent volume by 1.2.

Reagent	Volume (μl)
NuPCR Gene Expression Master Mix	10
NuPCR Gene Expression Assays Mix 1	1
NuPCR Gene Expression Assays Mix 2	1
PCR Grade Water	enough to bring the total volume to 20 μl per sample
Total Volume per Sample	20 - DNA volume

- 2 Transfer the total volume of reaction mix per sample to a new PCR tube or to each well of a new PCR plate.
- 3 Add the appropriate volume of DNA sample to the PCR tube or to each well of the PCR plate that contains the reaction mix. The total volume per tube or plate well should be 20 μl per sample.
- 4 Cap the PCR tube or seal the PCR plate with an adhesive seal.
- 5 Briefly centrifuge the tube or plate to remove air bubbles.
- 6 Place the PCR tube or plate in the pre-programmed real-time PCR instrument. Close the lid and start the amplification.

Data Analysis

NuPCR data analysis varies depending on the qPCR instrument being used. Reference the appropriate qPCR instrument user guide for data analysis instructions.

Notes

Technical Assistance

For technical assistance, contact your NuPCR distributor.

Table 5 Illumina General Contact Information

Illumina Website	http://www.illumina.com
-------------------------	---

MSDSs

Material safety data sheets (MSDSs) are available on the Illumina website at <http://www.illumina.com/msds>.

Product Documentation

You can obtain PDFs of additional product documentation from the Illumina website. Go to <http://www.illumina.com/support> and select a product. To download documentation, you will be asked to log in to MyIllumina. After you log in, you can view or save the PDF. To register for a MyIllumina account, please visit <https://my.illumina.com/Account/Register>.

Illumina

Headquartered in San Diego, California, U.S.A.

+1.800.809.ILMN (4566)

+1.858.202.4566 (outside North America)

techsupport@illumina.com

www.illumina.com