

coverage of traditionally challenging areas such as high AT/GC-rich regions, promoters, and homopolymeric regions.¹¹ To see a complete list of Illumina library preparation kits, visit support.illumina.com/sequencing/kits.html.

To advance the process even further, Illumina has combined the precision of digital microfluidics with its ease-of-use principles to create NeoPrep™ Library Prep System—a complete, fully automated library preparation instrument. Automation of library preparation will reduce opportunities for error, increase reproducibility, and reduce the amount of hands-on time required for a process that is often a bottleneck in the sequencing workflow. For more information on library prep automation developments, visit www.illumina.com/systems.html.

Multiplexing

In addition to the rise of data output per run, the sample throughput per run in NGS has also increased over time. Multiplexing allows large numbers of libraries to be pooled and sequenced simultaneously during a single sequencing run (Figure 5). With multiplexed libraries, unique index sequences are added to each DNA fragment during library preparation so that each read can be identified and sorted before final data analysis. With PE sequencing and multiplexing, NGS has dramatically reduced the time to data for multi-sample studies and enabled researchers to go from experiment to data faster and easier than ever before.




Figure 5: Library Multiplexing Overview.

- Two distinct libraries are attached to unique index sequences. Index sequences are attached during library preparation.
- Libraries are pooled together and loaded into the same flow cell lane.
- Libraries are sequenced together during a single instrument run. All sequences are exported to a single output file.
- A demultiplexing algorithm sorts the reads into different files according to their indexes.
- Each set of reads is aligned to the appropriate reference sequence.

Flexible, Scalable Instrumentation

While the latest NGS platforms can produce massive data output, NGS technology is also highly scalable. Sequencing systems are available for every method and scale of study, from small laboratories to large genome centers (Figure 6). Illumina NGS instruments range from the desktop MiSeq® Series, with output ranging from 0.3–15 Gb for small genome, amplicon, or targeted sequencing studies, to the colossal HiSeq X Ten fleet, which can generate an impressive, 16–18 Tb per run* for population-scale studies.

* With the full suite of 10 HiSeq X Systems.

