

Examples of Allelic Imbalance in Cell Lines


The ability of the HumanHap550 Genotyping BeadChip to detect allelic imbalances in the cancer cell lines HL60 (promyelocytic leukemia) and MCF7 (breast adenocarcinoma) is shown in Figure 2. All of the plots shown were generated with the Illumina Genome Viewer in BeadStudio. In HL60, deletions of varying size appear on Chromosome 9: a ~20Mb deletion on the p arm, and a ~1.8Mb deletion on the q arm (Figure 2A). These hemizygous deletions (from two copies to one copy) are manifest as a downward deflection in the log R ratio and a loss of heterozygotes in the AF. Other examples of aberrations in HL60 include monoallelic duplication of Chromosome 18 as indicated by an increase in the log R ratio and the split of the heterozygotes into two states: one at 0.67 (2:1 ratio) AF and another at 0.33 (1:2 ratio) AF (Figure 2B). There are also several small amplifications of approximately 260kb in HL60 on Chromosome 8 with 1:9 and 9:1 allelic ratios (Figure 2C). All of these aberrations have previously been identified using SKY karyotyping⁸. Adjacent homozygous and hemizygous deletions are revealed in MCF7 (Figure 2D). The gene annotation provided

in the Illumina BeadStudio software allows users to quickly analyze affected gene regions. In this case, both copies of MTAP are homozygously deleted and CDKN2A (p16) and CDKN2B (p15) are hemizygously deleted.

Summary

Illumina's Infinium Assays and Whole-Genome Genotyping BeadChips offer genome-wide coverage at an unprecedented resolution. The combined measurement of allelic ratios and normalized intensities provides enhanced detection of aberrations while facilitating identification of copy-neutral genetic anomalies such as uniparental disomy (UPD) and mitotic recombination^{9,10}. Whole-Genome Genotyping can also yield allelic information on deletions, duplications, and amplifications, which have implications in cancer therapeutics. The combination of proven assays, high-density arrays, and integrated software enables the analysis of LOH and copy number changes in both single and paired (i.e., matched) samples with high precision.

Figure 2: Detection of Allelic Imbalances in the Cancer Cell Lines HL60 and MCF7



- Genomic profiles of cancer cell lines on the HumanHap550 Genotyping BeadChip generated with Illumina BeadStudio software.
(A) Two regions in HL60 that exhibit LOH are shown by a downward deflection in the log R ratio and loss of heterozygotes in the AF profile.
(B) The entire length of Chromosome 18 in HL60 exhibits a monoallelic duplication (trisomy), as evidenced by the elevated log R ratio and split heterozygous cluster.
(C) Several amplifications in HL60 are present on Chromosome 8.
(D) High-resolution analysis of MCF7 showing a 150kb homozygous and 83kb hemizygous deletion in a region on Chromosome 9 containing MTAP, CDKN2A (p16), and CDKN2B (p15). For all log R ratio profiles, the blue line indicates a 100kb moving median for the HumanHap550 Genotyping BeadChip.

AAAGAATGATAACAGTAACACACCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAACGAACGATCAATTGAGACTAAATATAACGTACCATTAAAGAGTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATTCA

