

Amplify DNA (Pre-Amp)

- 1 Select **MSA8 Tasks | Make MSA8**.
 - a Select the WG#-DNA plate type.
- 2 Add reagents to quarter reservoirs (these volumes are for 3 plates only):

Reagent	Volume
MA1	9 ml
0.1 N NaOH	5 ml
MA2	13.5 ml
MSM	15 ml

- 3 Place the WG#-DNA plates and MSA8 plates on the robot deck.
- 4 Select **Run**.
- 5 Vortex the MSA8 plates at 1600 rpm for 1 minute.
- 6 Centrifuge the MSA8 plates at 280 × g at room temperature for 1 minute.

Incubate DNA

- 1 Select **Infinium HTS Extra | Incubate MSA8**.
 - a Scan the barcode of each MSA8 plate.
- 2 Incubate the MSA8 plates for 20–24 hours at 37°C.

Fragment DNA

- 1 Centrifuge the MSA8 plates at 280 × g at room temperature for 1 minute.
- 2 Select **MSA8 Tasks | Fragment MSA8**.
- 3 Place six MSA8 plates on the robot deck.
- 4 Add 20 ml FMS to a quarter reservoir.
- 5 Select **Run**.
- 6 **[Optional]** Adjust the Tecan scanner bracket to **Position A**.
- 7 Scan the barcode of the reagent bottle.
- 8 Select **OK**.
- 9 Vortex at 1600 rpm for 1 minute.
- 10 Centrifuge at 280 × g at room temperature for 1 minute.
- 11 Incubate at 37°C for 30 minutes.

SAFE STOPPING POINT

If you are stopping, seal the plates, and store at -25°C to -15°C.

Precipitate DNA

- 1 Select **MSA8 Tasks | Precip MSA8**.
- 2 Place six MSA8 plates on the robot deck.
- 3 Add 40 ml PM1 to a quarter reservoir.
- 4 Add 150 ml 2-propanol to a full reservoir.
- 5 Select **Run**.
 - a Scan the barcode of the reagent bottle.
 - b Select **OK**.
- 6 Invert the plates 10 times.
- 7 Select **Infinium HTS Extra | Spin MSA8**.
 - a Scan the barcode of each MSA8 plate.
- 8 Centrifuge at 3000 × g at 4°C for 20 minutes.
- 9 Invert the plates, and drain the supernatant.
- 10 Tap the plates several times.
- 11 Air dry for 15 minutes.

SAFE STOPPING POINT

If you are stopping, seal the plates, and store at -25°C to -15°C.

Resuspend DNA

- 1 Select **MSA8 Tasks | Resuspend MSA8**.
- 2 Place six MSA8 plates on the robot deck.
- 3 Add 20 ml RA1 to a quarter reservoir.
- 4 Select **Run**.
 - a Scan the barcode of the reagent bottle.
 - b Select **OK**.
- 5 Apply foil heat seals to the MSA8 plates.
- 6 Incubate for 15 minutes at 48°C.
- 7 Vortex at 1800 rpm for 1 minute.
- 8 Centrifuge at 280 × g for 1 minute.

SAFE STOPPING POINT

If you are stopping, store sealed MSA8 plates at 2°C to 8°C for up to 24 hours. If more than 24 hours, store at -25°C to -15°C. Store sealed RA1 at -25°C to -15°C. If RA1 will be used the next day, seal it, and store it overnight at 4°C.

Hybridize to BeadChip

- 1 Incubate the MSA8 plates at 95°C for 20 minutes.
- 2 Cool at room temperature for 30 minutes.
- 3 Centrifuge at 1500 × g at room temperature for 1 minute.
- 4 Place the gasket into the hyb chamber.
- 5 Dispense 400 µl PB2 into each reservoir.
- 6 Place the hyb chamber insert into the hyb chamber.
- 7 Close the hyb chamber.
- 8 Select **Infinium HTS Extra | Confirm for Hyb**.
- 9 Scan the barcodes.
- 10 Remove all BeadChips from packaging.
- 11 Place BeadChips into the robot BeadChip alignment fixtures.
- 12 Select **MSA8 Tasks | Hyb**.
 - a Select the 24-sample BeadChip.
 - b Enter the **Number of MSA8 plates**.
- 13 Place the robot BeadChip alignment fixtures onto the robot deck.
- 14 Place the MSA8 plates onto the robot deck.
- 15 Select **Run**.
- 16 Place each robot tip alignment guide on top of each robot BeadChip alignment fixture.
- 17 Select **OK**.
- 18 Select **OK**.
- 19 Remove the robot BeadChip alignment fixtures.
- 20 Place each BeadChip in a hyb chamber insert.
- 21 Close the hyb chamber.
- 22 Select **Infinium HTS Extra | Prepare Hyb Chamber**.
 - a Scan the barcodes.
- 23 Incubate at 48°C for 16 to 24 hours.

Prepare for Next Day

- 1 Soak the robot tip alignment guides in 1% aqueous Alconox solution.
- 2 Rinse and dry the robot tip alignment guides.
- 3 Add 330 ml 100% EtOH to the XC4 bottle and shake.

Wash BeadChips

- 1 Submerge the wash rack in the 1X PB1 wash.
- 2 Remove the hyb chamber inserts.
- 3 Inspect the BeadChips.
- 4 Remove BeadChips from the hyb chamber inserts.
- 5 Remove the cover seals from the BeadChips.
- 6 Place the BeadChips into the submerged wash rack.
- 7 Move the wash rack up and down for 1 minute.
- 8 Move the wash rack to the next 1X PB1 Wash.
- 9 Move the wash rack up and down for 1 minute.
- 10 Fill the multi-sample BeadChip alignment fixture with 150 ml 1X PB1.
- 11 Place black frames into the multi-sample BeadChip alignment fixture.
- 12 Place BeadChips into black frames.
- 13 Inspect the BeadChip. Remove excess residue.
- 14 Place a clear LCG spacer onto each BeadChip.
- 15 Place the alignment bar onto the multi-sample BeadChip alignment fixture.
- 16 Place LCG glass back plates on top of the clear spacers.
- 17 Attach the metal clamps to the flow-through chambers.
- 18 Trim the excess ends of the clear plastic spacers.
- 19 Return the flow-through chamber to the multi-sample BeadChip alignment fixture.
- 20 Select **Infinium HTS Extra | Wash**.
 - a Scan the barcode of the reagent bottles.

- b Scan the barcodes.

Extend and Stain (XStain)

- 1 Fill the water circulator.
- 2 Select **Robot QC Tasks | Circulator Manager** to set to 44°C.
- 3 Select **XStain Tasks | XStain LCG BeadChip HT**.
- 4 Add the following reagents to reservoirs:

Reagent	# BeadChips	Volume
95% formamide/1 mM EDTA	24	30 ml
	48	60 ml
RA1	24	30 ml
	48	60 ml
XC3	24	150 ml
	48	250 ml

- 5 Place the XStain plates on the robot deck.
- 6 Invert the EML tubes to mix, remove the caps, and place the EML tubes on the robot deck.
- 7 Enter the number of BeadChips.
 - a Select **Run**.
 - b Enter the stain temperature listed on the XStain plate.
 - c Select **OK**.
- 8 Place the LCG flow-through chambers into the chamber rack.
- 9 Select **OK**.
 - a Scan the barcode of the reagent bottle.
- 10 While the XStain task runs, wash the hybridization chamber humidifying buffer reservoirs.
- 11 Remove the LCG flow-through chambers from the chamber rack.
- 12 Set up PB1 and XC4 wash dishes.
- 13 Pour 310 ml PB1 into a wash dish.
- 14 Disassemble each LCG flow-through chamber.
- 15 Place BeadChips into a staining rack in the PB1 wash dish.
- 16 Submerge the LCG glass back plates in the DI H₂O wash basin.
- 17 Move the staining rack up and down 10 times.
- 18 Soak the BeadChips for 5 minutes.
- 19 Shake the XC4 bottle vigorously.
- 20 Pour 310 ml XC4 into a wash dish.
- 21 Move the staining rack to the XC4 wash dish.
- 22 Move the staining rack up and down 10 times.
- 23 Soak the BeadChips for 5 minutes.
- 24 Remove the staining rack.
- 25 Dry the BeadChips for 50–55 minutes at 675 mm Hg (0.9 bar).
- 26 Turn on the iScan™ systems.
- 27 Image the BeadChips immediately, or store them, protected from light.
- 28 To begin imaging, select **Infinium HTS Extra | Coat**.
 - a Scan the barcodes.