

For Research Use Only. Not for use in diagnostic procedures.

## Amplify DNA

| □ 1         | Add DNA into either of the following to create a             |
|-------------|--|
|             | DNA plate:   |
|             | Midi plate: 20 μl to each DNA well                           |
|             | ► TCY plate: 10 µl to each DNA well                          |
| $\square$ 2 | Select MSA3 Tasks   Make MSA3.                               |
| □3          | Select the DNA plate type.                                   |
| $\Box 4$    | Enter the <b>Number of DNA plates</b> .                      |
| □ 5         | •  |
|             | robot tube rack.   |
| □6          | Pour 15 ml NaOH into a trough and place on                   |
|             | the robot bed.   |
| $\square$ 7 | Place DNA and MSA3 plates on robot bed.                      |
| 8           | Select Run.  |
| 9           | Enter the barcode of each DNA plate.                         |
| □ 10        | Place the DNA plates on the robot bed and select <b>OK</b> . |
| □11         | Vortex the MSA3 plate at 1600 rpm for                        |
|             | 1 minute.  |
| □ 12        | Centrifuge at 280 × g at 22°C for 1 minute.                  |
| □ 13        | Remove the cap mat, place the MSA3 plate                     |
|             | on the robot bed, and select <b>OK</b> .                     |
| □ 14        | When complete, select OK.                                    |
| □ 15        | Remove and seal the MSA3 plate.                              |
| □16         | Centrifuge at 280 × g.                                       |

### Incubate DNA

| $\Box$ 1 | [LIMS] Select Infinium HTS                       |
|----------|--|
|          | a Scan the barcodes.                             |
| □ 2      | Incubate the MSA3 plate for 20–24 hours at 37°C. |

## Fragment DNA

| □ 1         | Centrifuge the MSA3 plate at 50 × g at room  |
|-------------|--|
|             | temperature for 1 minute.                    |
| $\square$ 2 | Select MSA3 Tasks   Fragment MSA3.           |
| □3          | Place the MSA3 plate on the robot bed.       |
| 4           | Place FMS tubes in the robot tube rack.      |
| $\Box$ 5    | Select Run.                                  |
| □6          | When complete, select OK.                    |
| □ 7         | Remove the plate and seal with a cap mat.    |
| □8          | Vortex at 1600 rpm for 1 minute.             |
| 9           | Centrifuge at 50 × g at room temperature for |
|             | 1 minute.                                    |
| □ 10        | Place on the 37°C heat block for 1 hour.     |

#### SAFE STOPPING POINT

If you are stopping, seal the plate(s), and store at -25°C to -15°C.



For Research Use Only. Not for use in diagnostic procedures.

### Precipitate DNA

| □ 1         | Select MSA3 Tasks   Precip MSA3.                |
|-------------|---|
| $\square$ 2 | Place the MSA3 plate on the robot bed.          |
| □3          | Place a half reservoir in the frame, and add    |
|             | 1 tube PM1 for 96 samples.                      |
| □ 4         | Place a full reservoir in the frame, and add 32 |
|             | ml 2-propanol.                                  |
| $\Box$ 5    | Select Run.                                     |
| □6          | Remove the MSA3 plate from the robot bed.       |
|             | Do not select <b>OK</b> .                       |
| $\Box$ 7    | Vortex at 1600 rpm for 1 minute.                |
| □8          | Incubate on the heat block for 5 minutes.       |
| 9           | Centrifuge at 50 × g for 1 minute.              |
| □10         | Set the centrifuge at 4°C.                      |
| □ 11        | Place the MSA3 plate on the robot bed.          |
| □ 12        | Select OK.                                      |
| □ 13        | Remove the MSA3 plate from the robot bed        |
|             | and seal.                                       |
| □ 14        | Invert 10 times to mix.                         |
| □ 15        | Incubate at 4°C for 30 minutes.                 |
| □16         | [LIMS] Select Infinium HTS   Spin MSA3.         |
| □ 17        | Select Run.                                     |
| □18         | Place in the centrifuge.                        |
| □ 19        | Centrifuge at 3000 × g for 20 minutes.          |
| □ 20        | Remove MSA3 plate.                              |
| □ 21        | Make sure that a blue pellet is present.        |
| 22          | Remove and discard the cap mat.                 |
| □ 23        | Quickly invert the plate and drain the          |
|             | supernatant.                                    |
| □ 24        | Firmly tap until all wells are free of liquid.  |
| □ 25        | Place the plate on a tube rack for 1 hour at    |
|             | room temperature.                               |
| □ 26        | Make sure that a blue pellet is still present.  |

#### SAFE STOPPING POINT

If you are stopping, seal the plate(s), and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$ .

### Resuspend DNA

| □ 1  | Select MSA3 Tasks   Resuspend MSA3.             |
|------|---|
| 2    | Place the MSA3 plate on the robot bed.          |
| □ 3  | Place a quarter reservoir in the frame, and     |
|      | add 7 ml RA1 for 96 samples.                    |
| 4    | Select Run.                                     |
| □ 5  | Remove the MSA3 plate from the robot deck.      |
| 6    | Apply a foil seal to the MSA3 plate.            |
| □ 7  | Incubate in the Illumina Hybridization Oven for |
|      | 1 hour.   |
| 8    | Vortex at 1800 rpm for 1 minute.                |
| 9    | Make sure that the pellets are resuspended.     |
| □ 10 | Pulse centrifuge at 280 × g.                    |

#### SAFE STOPPING POINT

If you are stopping, seal the plate, and store at -25°C to -15°C for up to 24 hours. For more than 24 hours, store at -80°C.

Store RA1 at -25°C to -15°C.



For Research Use Only. Not for use in diagnostic procedures.

| Hyb         | ridize DNA to the BeadChip   | ☐ 23 Incubate at 48°C for 16-24 hours |
|-------------|--|---------------------------------------|
| □1          | Incubate the MSA3 plate on the heat block for 20 minutes.                                    |                                       |
| $\square$ 2 | Cool at room temperature for 30 minutes.   |                                       |
| $\square$ 3 | Pulse centrifuge at 280 × g.   |                                       |
| □ 4         | Place the gasket into the hybridization chamber.   |                                       |
| $\square$ 5 | Add 400 µl PB2 into each reservoir.  |                                       |
| □ 6         | Place the hybridization chamber insert into the  |                                       |
| □ 7         | hybridization chamber.   |                                       |
| ∐7<br>□8    | Immediately cover the chamber with the lid.  [LIMS] Select Select Infinium HTS   Confirm for |                                       |
|             | Hyb.   |                                       |
| 9           | [LIMS] Scan the barcodes.  |                                       |
| □ 10        | Remove all BeadChips from packaging.   |                                       |
| □ 11        | Place BeadChips into the robot BeadChip  |                                       |
|             | alignment fixtures.  |                                       |
| □12         | Place the robot BeadChip alignment fixtures onto the robot deck.                             |                                       |
| □13         | Pulse centrifuge the MSA3 plate at 280 × g.  |                                       |
|             | Place the MSA3 plate onto the robot deck.  |                                       |
|             | Select Run.  |                                       |
|             | Place each robot tip alignment guide on top of   |                                       |
|             | each robot BeadChip alignment fixture.   |                                       |
| □ 17        | To start the run, select <b>OK</b> .   |                                       |
|             | When complete, select <b>OK</b> .  |                                       |
| □19         | Remove the robot BeadChip alignment  |                                       |
|             | fixtures.  Place each BeadChip in a hybridization  |                                       |
| L 20        | chamber insert.  |                                       |
| □21         | Place the lid on the chamber and secure with   |                                       |
|             | the metal clamps.  |                                       |
| □ 22        | [LIMS] Select Infinium HTS   Prepare Hyb   |                                       |
| _           | Chamber.   |                                       |
|             | a Scan the barcodes.   |                                       |

## Prepare for Next Day

| □ 1         | Add 330 ml fresh 100% EtOH to the XC4         |
|-------------|---|
|             | bottle.                                       |
| $\square$ 2 | Vigorously shake to resuspend.                |
| □ 3         | Leave the bottle upright on the lab bench     |
|             | overnight.                                    |
| □ 4         | Soak the robot tip alignment guides in 1%     |
|             | aqueous Alconox solution.                     |
| □ 5         | Rinse and dry the robot tip alignment guides. |



For Research Use Only. Not for use in diagnostic procedures.

| Wash BeadChips       |   | with DI H <sub>2</sub> O. |          | Extend and Stain BeadChips   |                   |             |  |
|----------------------|---|---------------------------|----------|--|-------------------|-------------|--|
| □1<br>□2<br>□3<br>□4 | Submerge the wash rack in the PB1 wash. Remove the hybridization insert. Remove the BeadChips. Remove the cover seals from the BeadChips. | WILLIAM TO THE CO.        |          | Fill the water circul<br>Select Robot QC T<br>to set to 44°C.<br>Select XStain Task                | asks   Circulator |             |  |
| 5                    | ·   |                           | □3<br>□4 | Add the following r  |                   |             |  |
|                      | Place the BeadChips into the submerged wash rack.   |                           |          | Reagent  | # BeadChips       | Volume      |  |
| ☐ 6                  | Move the wash rack up and down for 1 minute.  |                           |          | 95% formamide/1<br>mM EDTA   | 1–8               | 15 ml       |  |
| 7                    | Move the wash rack to the next PB1 Wash.  |                           |          |  | 9–16              | 17 ml       |  |
| 8                    | Move the wash rack up and down for 1  |                           |          |  | 17–24             | 25 ml       |  |
| 9                    | minute.  Confirm that you are using the correct Infinium  |                           |          |  | 1–8               | 10 ml       |  |
|                      | LCG glass back plates and spacers.  |                           |          |  | 9–16              | 20 ml       |  |
| <u> </u>             | Fill the BeadChip alignment fixture with 150 ml   |                           |          |  | 17–24             | 30 ml       |  |
|                      | PB1.  |                           |          | \\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\   |                   |             |  |
| 11                   | For each BeadChip, place one black frame  |                           |          | XC3  | 1–8               | 50 ml       |  |
|                      | into the BeadChip alignment fixture.  |                           |          |  | 9–16              | 100 ml      |  |
|                      | Place each BeadChip into a black frame.   |                           |          |  | 17–24             | 150 ml      |  |
| ∟13                  | Place a <i>clear</i> spacer onto the top of each  |                           | □5       | Invert the LX1 LX2   | P FML SML an      | d ATM       |  |
| □ 14                 | BeadChip. Place the alignment bar onto the alignment fixture.   |                           |          | Invert the LX1, LX2, EML, SML, and ATM tubes to mix. Remove the caps, and place on the robot deck. |                   |             |  |
| 15                   | Place a clean glass back plate on top of each   |                           | □6       | Enter the number of  | of BeadChips.     |             |  |
|                      | clear spacer.   |                           | $\Box$ 7 | Select Run.  |                   |             |  |
| □16                  | Secure each flow-through chamber assembly with metal clamps.  |                           | □8       | [Non-LIMS]Enter the stain temperature listed on the SML tube.                                      |                   |             |  |
| <u> </u>             | Remove the assembled flow-through chamber   |                           | □9       | Place the flow-thro  | ugh chambers in   | nto the     |  |
|                      | from the alignment fixture.   |                           |          | chamber rack.  |                   |             |  |
| □ 18                 | Trim the spacers from each end of the   |                           |          | Select <b>OK</b> .   |                   |             |  |
|                      | assembly.   |                           | □11      | Remove the flow-ti   | hrough chambei    | rs from the |  |
| <u> </u>             | Leave assembled flow-through chambers on  |                           |          | chamber rack.  |                   |             |  |
| the lab bench.       |   |                           | ∟12      | Set up two top-loa PB1 and XC4.  | ding wash dishe   | s labeled   |  |



For Research Use Only. Not for use in diagnostic procedures.

| □13          | Add 310 ml PB1 to the PB1 wash dish.                        |
|--------------|---|
| □ 14         | Submerge the staining rack in the wash dish.                |
| □ 15         | Leave the staining rack in the wash dish.                   |
| □16          | Disassemble each flow-through chamber.                      |
| □ 17         | Place the BeadChips into the submerged                      |
|              | staining rack.  |
| □18          | Slowly lift the staining rack 10 times.                     |
| □ 19         | Soak for 5 minutes.   |
| $\square$ 20 | Vigorously shake the XC4 bottle.                            |
| □ 21         | Add 310 ml XC4 to the XC4 wash dish and                     |
|              | cover.  |
| $\square$ 22 | Transfer the staining rack from the PB1 to the              |
|              | XC4.  |
|              | Slowly lift the staining rack 10 times.                     |
| _            | Soak for 5 minutes.   |
| □ 25         | Remove the staining rack and place it onto the              |
| _            | tube rack.  |
|              | Dry each BeadChip as follows.                               |
|              | a Grip the BeadChip by the barcode end.                     |
|              | b Place onto a tube rack with the barcode                   |
| □ 07         | facing up and toward you.                                   |
| L 21         | Place the tube rack into the vacuum desiccator.             |
| □ 28         |   |
| ∟20          | Dry the BeadChips for 50–55 minutes at 675 mm Hg (0.9 bar). |
| □ 20         | [LIMS] Select Infinium HTS   Coat.                          |
|              | a Scan the barcodes.  |
|              | a countilo balocaco.  |