# illumina

## Infinium HD Methylation Assay Manual Workflow Checklist

## Convert DNA

- □ 1 Follow the instructions in the Zymo EZ DNA Methylation Kit to denature the genomic DNA and add conversion reagent.
- 2 Incubate in a thermal cycler using the following settings for 16 cycles:
  - ▶ 95°C for 30 seconds
  - ▶ 50°C for 1 hour
- □ 3 Hold DNA at 4°C for 10 minutes until cleanup.
- ☐ 4 Use the instructions in the Zymo EZ DNA Methylation Kit to cleanup the conversion reagent.

#### SAFE STOPPING POINT

If you are stopping, heat-seal the plate, and store at -25°C to -15°C for up to 30 days.

# Create the BCD Plate

- □ 1 If frozen, thaw BCD samples to room temperature and vortex to mix.
- 2 Apply a BCD barcode label to a new 0.8 ml midi plate or a new 0.2 ml TCY plate.
- □ 3 Transfer the BCD to the plate as follows
  - Midi plate: 20 µl BCD sample to each well (requires ≥ 1000 ng input in bisulfite conversion)
  - TCY plate: 10 µl BCD sample to each well (requires ≥ 500 ng input in bisulfite conversion)

## Amplify DNA

- □ 1 Add 20 µl MA1 to each well.
- 2 Transfer 4 µl DNA sample from the BCD plate to the MSA4 plate.
- $\Box$  3 Add 4 µl 0.1N NaOH in to each well.
- $\Box$  4 Seal the MSA4 plate with the 96-well cap mat.
- $\Box$  5 Vortex at 1600 rpm for 1 minute, and then pulse centrifuge at 280 × g.
- 6 Incubate at room temperature for 10 minutes.
- $\Box$  7 Add 68 µl RPM in to each well.
- $\square$  8 Add 75 µl MSM in to each well.
- □ 9 Vortex at 1600 rpm for 1 minute, and then pulse centrifuge at 280 × g.

Incubate DNA

□ 1 Incubate the MSA4 plate for 20–24 hours at 37°C.

## Fragment DNA

- $\Box$  1 Pulse centrifuge the plate at 280 × g.
- $\Box$  2 Add 50 µl FMS to the MSA4 plate.
- □ 3 Vortex at 1600 rpm for 1 minute, and then centrifuge at  $280 \times g$  for 1 minute.
- $\Box$  4 Incubate at 37° C for 1 hour.

#### SAFE STOPPING POINT

If you are stopping, seal the plate, and store at -25°C to -15°C for up to 24 hours.

## Precipitate DNA

- □ 1 Add 100 µl PM1 to the MSA4 plate.
- $\Box$  2 Reseal with the cap mat.
- □ 3 Vortex the plate at 1600 rpm for 1 minute.
- $\Box$  4 Incubate at 37° C for 5 minutes.
- $\Box$  5 Pulse centrifuge at 280 × g for 1 minute.
- $\Box$  6 Set the centrifuge at 4°C.
- $\Box$  7 Remove and discard the cap mat.
- $\square\,8$  Add 300  $\mu l$  100% 2-propanol to each well.
- $\Box$  9 Apply fresh cap mats.
- $\Box$  10 Invert the plate 10 times to mix.
- □ 11 Incubate in a refrigerator set at 4°C for 30 minutes.
- $\Box$  12 Centrifuge at 3000 × g at 4°C for 20 minutes.
- 13 Immediately remove the plate from the centrifuge.
- $\Box$  14 Make sure that a blue pellet is present.
- $\Box$  15 Remove and discard the cap mat.
- 16 Quickly invert the plate and drain the supernatant.
- $\Box$  17 Firmly tap until all wells are free of liquid.
- □ 18 Place the plate on the tube rack for 1 hour at room temperature.
- $\Box$  19 Make sure that a blue pellet is still present.

### SAFE STOPPING POINT

If you are stopping, seal the plate, and store at -25°C to -15°C for up to 24 hours.

## **Resuspend DNA**

- $\Box$  1 Add 46 µl RA1 per well.
- $\Box$  2 Apply a foil heat seal.
- □ 3 Incubate at 48°C for 1 hour.
- $\Box$  4 Vortex at 1800 rpm for 1 minute.
- $\Box$  5 Pulse centrifuge at 280 × g.

#### SAFE STOPPING POINT

If you are stopping, seal the plate, and store at -25°C to -15°C for up to 24 hours.

- Hybridize DNA to the BeadChip
- □ 1 Incubate the MSA4 plate at 95° C on the heat block for 20 minutes.
- $\Box$  2 Cool at room temperature for 30 minutes.
- $\Box$  3 Pulse centrifuge at 280 × g.
- 4 Place the gasket into the hybridization chamber.
- $\Box$  5 Add 400 µl PB2 to the top and bottom wells.
- $\Box$  6 Immediately cover the chamber with the lid.
- $\Box$  7 Remove the BeadChips from packaging.
- $\square$  8 Place each BeadChip into an insert.
- $\Box$  9 Transfer 26 µl each sample to the BeadChip.
- $\Box$  10 Wait for the DNA to disperse.
- $\Box$  11 Inspect the loading port for excess liquid.
- 12 If excess liquid is not present, add leftover sample.
- □ 13 Store RA1 at -25°C to -15°C.
- 14 Load the inserts into the hybridization chamber.
- □ 15 Place the lid on the chamber and secure with the metal clamps.
- ☐ 16 Place the chamber into the preheated Illumina Hybridization Oven.
- $\Box$  17 Incubate at 48°C for 16–24 hours.
- $\Box$  18 Store RA1 at -25°C to -15°C.

## Prepare for Next Day

- □ 1 Add 330 ml fresh 100% EtOH to the XC4 bottle.
- 2 Leave the bottle upright on the lab bench overnight.

# illumina

## Infinium HD Methylation Assay Manual Workflow Checklist

## Wash BeadChips

- $\Box$  1 Submerge the wash rack in the PB1 wash.
- $\Box$  2 Remove the hybridization insert.
- $\Box$  3 Remove the BeadChips.
- $\Box$  4 Remove the cover seals from the BeadChips.
- □ 5 Place the BeadChips into the submerged wash rack.
- 6 Move the wash rack up and down for 1 minute.
- $\Box$  7 Move the wash rack to the next PB1 Wash.
- □8 Move the wash rack up and down for 1 minute.
- 9 Confirm that you are using the correct Infinium glass back plates and spacers.
- □ 10 Fill the BeadChip alignment fixture with 150 ml PB1.
- 11 For each BeadChip, place one black frame into the BeadChip alignment fixture.
- $\Box$  12 Place each BeadChip into a black frame.
- □ 13 Place a *clear* spacer onto the top of each BeadChip.
- 14 Place the alignment bar onto the alignment fixture.
- 15 Place a clean glass back plate on top of each clear spacer.
- ☐ 16 Secure each flow-through chamber assembly with metal clamps.
- 17 Remove the assembled flow-through chamber from the alignment fixture.
- 18 Trim the spacers from each end of the assembly.
- 19 Leave assembled flow-through chambers on the lab bench.

20 Wash the hybridization chamber reservoirs with DI H<sub>2</sub>O.

## Extend and Stain BeadChips

- $\Box$  1 Fill the water circulator.
- 2 Turn on the water circulator and set the temperature to 44°C.
- □ 3 When the chamber rack reaches 44°C, place the flow-through chamber assemblies into the chamber rack.
- 4 Fill the reservoir of each flow-through chamber as follows.
  - a 150 µl RA1. Incubate for 30 seconds. Repeat 4 times.
    - [\_] 1 [\_] 2 [\_] 3 [\_] 4 [\_] 5
  - $\Box$  b 450 µl XC1. Incubate for 10 minutes.
  - □ c 450 µl XC2. Incubate for 10 minutes.
  - d 200 µl TEM. Incubate for 15 minutes.
  - e 450 µl 95% formamide/1 mM EDTA. Incubate for 1 minute. Repeat once.
    - [\_] 1 [\_] 2
  - $\Box$  f Incubate 5 minutes.
  - g Set the the chamber rack temperature to the temperature indicated on the STM tube.
  - h 450 µl XC3. Incubate for 1 minute. Repeat once.
    - [\_] 1 [\_] 2
- □ 5 Wait for the chamber rack to reach the correct temperature.
- 6 If imaging the BeadChip immediately after the staining process, turn on the scanner.
- ☐ 7 Fill the reservoir of each flow-through chamber as follows.
  - $\Box$  a 250 µl STM. Incubate for 10 minutes.
  - b 450 µl XC3. Incubate for 1 minute. Repeat once.

# illumina

## Infinium HD Methylation Assay Manual Workflow Checklist

[\_] 1 [\_] 2

- C Wait 5 minutes.
- d 250 µl ATM. Incubate for 10 minutes.
- e 450 µl XC3. Incubate for 1 minute. Repeat once.

### []1[]2

- $\Box$  f Wait 5 minutes.
- g 250 µl STM. Incubate for 10 minutes.
- 450 µl XC3. Incubate for 1 minute. Repeat once.

### [\_] 1 [\_] 2

- □ i Wait 5 minutes.
- i 250 µl ATM. Incubate for 10 minutes.
- k 450 µl XC3. Incubate for 1 minute. Repeat once.

### [\_] 1 [\_] 2

- □ I Wait 5 minutes.
- m 250 µl STM. Incubate for 10 minutes.
- n 450 µl XC3. Incubate for 1 minute. Repeat once.

## [\_] 1 [\_] 2

- $\Box$  o Wait 5 minutes.
- 8 Remove the flow-through chambers from the chamber rack.
- 9 Set up two top-loading wash dishes labeled PB1 and XC4.
- 10 Add 310 ml PB1 to the PB1 wash dish.
- $\Box$  11 Submerge the staining rack in the wash dish.
- $\Box$  12 Leave the staining rack in the wash dish.
- □ 13 Disassemble each flow-through chamber.
- 14 Place the BeadChips into the submerged staining rack.
- $\Box$  15 Slowly move the staining rack up and down 10 times.

- $\Box$  17 Vigorously shake the XC4 bottle.
- 18 Add 310 ml XC4 to the XC4 wash dish and cover.
- 19 Transfer the staining rack to the XC4 wash dish.
- 20 Slowly lift the staining rack up and down 10 times.
- $\square$  21 Soak for 5 minutes.
- $\square$  22 Remove the staining rack and place it onto the tube rack.
- 23 Place the tube rack into the vacuum desiccator.
- 24 Dry the BeadChips for 50-55 minutes at 675 mm Hg (0.9 bar).

### SAFE STOPPING POINT

Store the BeadChips in the Illumina BeadChip Slide Storage Box at room temperature. Scan within 72 hours.