


Nextera DNA Library Prep Protocol Guide

For Research Use Only. Not for use in diagnostic procedures.

Tagment Genomic DNA	3
Clean Up Tagmented DNA	4
Amplify Tagmented DNA	5
Clean Up Libraries	6
Check Libraries	7
Normalize and Pool Libraries	8
Acronyms	9
Technical Assistance	11



ILLUMINA PROPRIETARY

Part # 1000000006836 v00
January 2016

illumina®

This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2016 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

Tagment Genomic DNA

Preparation

- 1 Label a new 96-well TCY plate NSP1.
- 2 Save the following program as TAG NSP1 on a thermal cycler with a heated lid.
 - ▶ Choose the preheat lid option and set to 55°C
 - ▶ 55°C for 5 minutes
 - ▶ Hold at 10°C

Procedure

- 1 Add 20 μl of genomic DNA at 2.5 $\text{ng}/\mu\text{l}$ (50 ng total) to each sample well of the NSP1 plate.
- 2 Add 25 μl of TD Buffer to the wells containing genomic DNA.
- 3 Add 5 μl of TDE1 to the wells containing genomic DNA and TD Buffer.
- 4 Pipette up and down 10 times to mix.
- 5 Centrifuge at $280 \times g$ at 20°C for 1 minute.
- 6 Place on the preprogrammed thermal cycler and run the TAG NSP1 program.

Clean Up Tagmented DNA

Preparation

- 1 Label a new midi plate NSP2.
- 2 Label a new TCY plate NSP3.
- 3 Add 180 μ l Zymo DNA binding buffer to each well of the NSP2 plate with a sample in the corresponding well of the NSP1 plate.

Procedure

- 1 Transfer 50 μ l from each well of NSP1 to the corresponding well of NSP2. Pipette up and down 10 times to mix.
- 2 Place the Zymo-Spin I-96 Plate on the Collection Plate.
- 3 Transfer sample mixture from NSP2 to the Zymo-Spin I-96 Plate.
- 4 Centrifuge at 1300 \times g at 20°C for 2 minutes.
- 5 Discard the flow-through.
- 6 Wash 2 times with 300 μ l Zymo wash buffer.
- 7 Centrifuge at 1300 \times g for 2 minutes.
- 8 Place the Zymo-Spin I-96 Plate on NSP3.
- 9 Add 25 μ l of RSB directly to the column matrix in each well.
- 10 Incubate for 2 minutes at room temperature.
- 11 Centrifuge at 1300 \times g at 20°C for 2 minutes.

Amplify Tagmented DNA

Preparation

- 1 Label a new 96-well microplate NAP1.
- 2 Save the following program as PCR AMP on a thermal cycler with a heated lid.
 - ▶ Choose the preheat lid option and set to 100°C
 - ▶ 72°C for 3 minutes
 - ▶ 98°C for 30 seconds
 - ▶ 5 cycles of:
 - ▶ 98°C for 10 seconds
 - ▶ 63°C for 30 seconds
 - ▶ 72°C for 3 minutes
 - ▶ Hold at 10°C

Procedure

- 1 Arrange Index 1 (i7) adapters in the TruSeq Index Plate Fixture as follows:
 - ▶ 24 libraries—Columns 1–6
 - ▶ 96 libraries—Columns 1–12
- 2 Arrange Index 2 (i5) adapters in the TruSeq Index Plate Fixture as follows:
 - ▶ 24 libraries—Rows A–D
 - ▶ 96 libraries—Rows A–H
- 3 Place the plate on the TruSeq Index Plate Fixture.
- 4 Using a multichannel pipette, add 5 μ l of each Index 1 (i7) adapter down each column. Replace the cap on each i7 adapter tube with a new orange cap.
- 5 Using a multichannel pipette, add 5 μ l of each Index 2 (i5) adapter across each row. Replace the cap on each i5 adapter tube with a new white cap.
- 6 Add 15 μ l NPM.
- 7 Add 5 μ l PPC.
- 8 Transfer 20 μ l from each well of NSP3 to NAP1. Pipette up and down 3–5 times to mix.
- 9 Centrifuge at $280 \times g$ at 20°C for 1 minute.
- 10 Transfer the NAP1 plate to the post-amplification area.
- 11 Place on the preprogrammed thermal cycler and run the PCR AMP program.

Clean Up Libraries

Preparation

- 1 Prepare fresh 80% ethanol from absolute ethanol.
- 2 Label a new midi plate NAP2.
- 3 Label a new TCY plate NLP.

Procedure

- 1 Centrifuge NAP1 at $280 \times g$ at 20°C for 1 minute.
- 2 Transfer the contents of NAP1 to NAP2.
- 3 Vortex AMPure XP beads for 30 seconds. Add beads to a trough.
- 4 Add $30 \mu\text{l}$ AMPure XP beads to NAP2.
- 5 Add $30 \mu\text{l}$ AMPure XP beads to NAP2.
- 6 Pipette up and down 10 times to mix.
- 7 Incubate at room temperature for 5 minutes.
- 8 Place on a magnetic stand and wait until the liquid is clear (~2 minutes). Keep on magnetic stand until step 13.
- 9 Remove and discard supernatant.
- 10 Wash 2 times with $200 \mu\text{l}$ 80% EtOH.
- 11 Remove residual EtOH.
- 12 Air-dry the beads for 15 minutes.
- 13 Remove from the magnetic stand.
- 14 Add $32.5 \mu\text{l}$ RSB.
- 15 Pipette up and down 10 times to mix.
- 16 Incubate at room temperature for 2 minutes.
- 17 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 18 Transfer $30 \mu\text{l}$ supernatant from NAP2 to NLP.

Check Libraries

Quantify Libraries

- 1 Quantify your libraries. Convert library concentration using the formula $1 \text{ ng}/\mu\text{l} = 3 \text{ nM}$.

Quality Control

- 1 Run $1 \mu\text{l}$ of 1:3 diluted library on an Agilent Technology 2100 Bioanalyzer using a High Sensitivity DNA chip.

Library Size from Bioanalyzer in bp	Conversion Factor for $\text{ng}/\mu\text{l} > \text{nM}$	DNA Concentration for Cluster Generation
250	$1 \text{ ng}/\mu\text{l} = 6 \text{ nM}$	6–12 pM
500	$1 \text{ ng}/\mu\text{l} = 3 \text{ nM}$	6–12 pM
1000–1500	$1 \text{ ng}/\mu\text{l} = 1.5 \text{ nM}$	12–20 pM

Normalize and Pool Libraries

Preparation

- 1 Apply the NDP barcode label to a new 96-well midi plate.
- 2 Apply the NPP barcode label to a new 96-well midi plate (for indexed libraries).
- 3 If NLP was stored frozen, thaw at room temperature and then centrifuge at $280 \times g$ for 1 minute.

Make NDP

- 1 Transfer 10 μl library from each well of NLP to the corresponding well of NDP.
- 2 Normalize the concentration in each well to 2 nM using Tris-Cl 10 mM, pH 8.5 with 0.1% Tween 20.
- 3 Shake at 1000 rpm for 2 minutes.
- 4 Centrifuge at $280 \times g$ for 1 minute.

Make NPP [For Indexed Libraries]

- 1 Transfer 5 μl from each well in column 1 of NDP to column 1 of NPP.
- 2 Repeat step 1 for the remaining columns of NDP until samples are pooled in column 1 of NPP.
- 3 Combine the contents of column 1 into A2 of NPP.
- 4 Shake at 1800 rpm for 2 minutes.
- 5 Denature and dilute pooled libraries to the loading concentration for the instrument you are using. See the denature and dilute libraries guide for your instrument.

Acronyms

Acronym	Definition
NAP1	Nextera Amplification Plate 1
NAP2	Nextera Amplification Plate 2
NDP	Nextera Dilution Plate
NLP	Nextera Library Plate
NPM	Nextera PCR Master Mix
NSP1	Nextera Sample Plate 1
NSP2	Nextera Sample Plate 2
NSP3	Nextera Sample Plate 3
NPP	Nextera Pooled Plate
PPC	PCR Primer Cocktail
RSB	Resuspension Buffer
TD	Tagment DNA Buffer
TDE1	Tagment DNA Enzyme

Notes

Technical Assistance

For technical assistance, contact Illumina Technical Support.

Table 1 Illumina General Contact Information

Website	www.illumina.com
Email	techsupport@illumina.com

Table 2 Illumina Customer Support Telephone Numbers

Region	Contact Number	Region	Contact Number
North America	1.800.809.4566	Japan	0800.111.5011
Australia	1.800.775.688	Netherlands	0800.0223859
Austria	0800.296575	New Zealand	0800.451.650
Belgium	0800.81102	Norway	800.16836
China	400.635.9898	Singapore	1.800.579.2745
Denmark	80882346	Spain	900.812168
Finland	0800.918363	Sweden	020790181
France	0800.911850	Switzerland	0800.563118
Germany	0800.180.8994	Taiwan	00806651752
Hong Kong	800960230	United Kingdom	0800.917.0041
Ireland	1.800.812949	Other countries	+44.1799.534000
Italy	800.874909		

Safety data sheets (SDSs)—Available on the Illumina website at support.illumina.com/sds.html.

Product documentation—Available for download in PDF from the Illumina website. Go to support.illumina.com, select a product, then select **Documentation & Literature**.



Illumina

5200 Illumina Way

San Diego, California 92122 U.S.A.

+1.800.809.ILMN (4566)

+1.858.202.4566 (outside North America)

techsupport@illumina.com

www.illumina.com