

## Tagment Genomic DNA

- 1 Quantify gDNA using a fluorometric method.
- 2 Dilute gDNA in Tris-HCl 10 mM, pH 8.5 to a final volume of 10  $\mu$ l at 5 ng/ $\mu$ l.
- 3 Add the following to a new plate or to a new tube.
  - ▶ TD (25  $\mu$ l)
  - ▶ Normalized gDNA (10  $\mu$ l)
  - ▶ TDE2 (15  $\mu$ l)
- 4 Mix thoroughly.
- 5 Centrifuge.
- 6 Place on the preprogrammed thermal cycler and run the TAG58 program.
- 7 Add 15  $\mu$ l ST2, and then pipette to mix
- 8 Place on the preprogrammed thermal cycler and run the TAG60 program.

## Clean Up Tagmented DNA

- 1 Transfer all supernatant.
- 2 Add 52  $\mu$ l SPB, and then mix thoroughly.
- 3 Incubate at room temperature for 5 minutes.
- 4 Place on a magnetic stand until the liquid is clear.
- 5 Transfer 98  $\mu$ l supernatant.
- 6 Add 137  $\mu$ l SPB, and then mix thoroughly.
- 7 Incubate at room temperature for 5 minutes.
- 8 Place on a magnetic stand until the liquid is clear.
- 9 Remove and discard all supernatant.
- 10 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 11 Using a 20  $\mu$ l pipette, remove residual 80% EtOH.
- 12 Air-dry on the magnetic stand for 5 minutes.
- 13 Add 22.5  $\mu$ l RSB, and then mix thoroughly.
- 14 Remove from the magnetic stand.
- 15 Incubate at room temperature for 2 minutes.
- 16 Centrifuge.
- 17 Place on a magnetic stand until the liquid is clear.
- 18 Transfer 20  $\mu$ l supernatant.

## Amplify Tagmented DNA

- 1 [Plate] Arrange Index 1 (i7) adapters in columns 1–12.
- 2 [Plate] Arrange Index 2 (i5) adapters in rows A–H.
- 3 [Plate] Place the plate on the TruSeq Index Plate Fixture.
- 4 Add 5  $\mu$ l of each Index 1 (i7) adapter.
- 5 Add 5  $\mu$ l of each Index 2 (i5) adapter.
- 6 Add 20  $\mu$ l LAM, and then mix thoroughly.
- 7 Centrifuge.
- 8 Place on the thermal cycler and run the LAM AMP program.

### SAFE STOPPING POINT

If you are stopping, seal the plate or cap the tube and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

## Clean Up Amplified DNA

- 1 Centrifuge.
- 2 Transfer 50  $\mu$ l total volume.
- 3 Add 90  $\mu$ l SPB, and then mix thoroughly.
- 4 Incubate at room temperature for 5 minutes.
- 5 Centrifuge.
- 6 Place on a magnetic stand until liquid is clear.
- 7 Remove and discard all supernatant.
- 8 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 9 Using a 20  $\mu$ l pipette, remove residual 80% EtOH.
- 10 Air-dry on the magnetic stand for 5 minutes.
- 11 Add 17  $\mu$ l RSB, and then mix thoroughly.
- 12 Remove from the magnetic stand.
- 13 Incubate at room temperature for 2 minutes.
- 14 Centrifuge.
- 15 Place on a magnetic stand until liquid is clear.
- 16 Transfer 15  $\mu$ l supernatant.
- 17 Quantify the library using a fluorometric method.

**SAFE STOPPING POINT**

If you are stopping, seal the plate or cap the tube and store at -25°C to -15°C for up to 14 days.

## Hybridize Probes

- 1 Combine 500 ng of each DNA library, making sure that each library has a unique index.
  - ▶ If the total volume is > 30  $\mu$ l, concentrate the pooled sample to 30  $\mu$ l.
  - ▶ If the total volume is < 30  $\mu$ l, increase the volume to 30  $\mu$ l with RSB.
- 2 Add the following to a new plate or to a new tube.
  - ▶ DNA library sample or pool (30  $\mu$ l)
  - ▶ BLR (10  $\mu$ l)
  - ▶ CEX (10  $\mu$ l)
- 3 Mix thoroughly.
- 4 Centrifuge.
- 5 Add 125  $\mu$ l SPB, and then mix thoroughly.
- 6 Incubate at room temperature for 10 minutes.
- 7 Centrifuge.
- 8 Place on a magnetic stand until the liquid is clear.
- 9 Remove and discard all supernatant.
- 10 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 11 Using a 20  $\mu$ l pipette, remove residual 80% EtOH.
- 12 Air-dry on the magnetic stand for 10 minutes.
- 13 Add 7.7  $\mu$ l EHB1, and then mix thoroughly.
- 14 Remove from the magnetic stand.
- 15 Incubate at room temperature for 2 minutes.
- 16 Centrifuge.
- 17 Place on a magnetic stand until the liquid is clear.
- 18 Transfer 7.5  $\mu$ l supernatant .
- 19 Add 2.5  $\mu$ l EHB2, and then mix thoroughly.
- 20 Centrifuge.
- 21 Place on the thermal cycler and run the TRE HYB program.

## Capture Hybridized Probes

- 1 Centrifuge.
- 2 Transfer all (~10  $\mu$ l).
- 3 Add 250  $\mu$ l SMB and mix.
- 4 Incubate at room temperature for 25 minutes.
- 5 Centrifuge.
- 6 Place on a magnetic stand until liquid is clear.
- 7 Remove and discard all supernatant.
- 8 Remove from the magnetic stand.
- 9 Add 200  $\mu$ l EEW and mix.
- 10 Incubate at 50°C as follows.
  - ▶ [Plate] Place on the microheating system for 30 minutes.
  - ▶ [Tube] Place on the heat block for 30 minutes.
- 11 Place on a magnetic stand until liquid is clear.
- 12 Remove and discard all supernatant.
- 13 Remove from the magnetic stand.
- 14 Repeat steps 9–13 for a total of 2 washes.
- 15 Mix 28.5  $\mu$ l EE1 and 1.5  $\mu$ l HP3, and then vortex.
- 16 Add 23  $\mu$ l elution premix and mix.
- 17 Incubate at room temperature for 2 minutes.
- 18 Centrifuge.
- 19 Place on a magnetic stand until liquid is clear.
- 20 Transfer 21  $\mu$ l supernatant.
- 21 Add 4  $\mu$ l ET2 and mix.
- 22 Add 5  $\mu$ l RSB and mix.
- 23 Centrifuge.

**SAFE STOPPING POINT**

If you are stopping, seal the plate or cap the tube and store at -25°C to -15°C for up to 7 days.

## Perform Second Hybridization

- 1 Add the following.
  - ▶ BLR (10 µl)
  - ▶ CEX (10 µl)
- 2 Mix thoroughly.
- 3 Centrifuge.
- 4 Add 125 µl SPB, and then mix thoroughly.
- 5 Incubate at room temperature for 10 minutes.
- 6 Centrifuge.
- 7 Place on a magnetic stand until the liquid is clear.
- 8 Remove and discard all supernatant.
- 9 Wash 2 times with 200 µl 80% EtOH.
- 10 Using a 20 µl pipette, remove residual 80% EtOH.
- 11 Air-dry on the magnetic stand for 10 minutes.
- 12 Add 7.7 µl EHB1, and then mix thoroughly.
- 13 Remove from the magnetic stand.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge.
- 16 Place on a magnetic stand until the liquid is clear.
- 17 Transfer 7.5 µl supernatant .
- 18 Add 2.5 µl EHB2, and then mix thoroughly.
- 19 Centrifuge.
- 20 Place on the thermal cycler and run the TRE HYB program.

## Perform Second Capture

- 1 Centrifuge.
- 2 Transfer 10 µl supernatant.
- 3 Add 250 µl SMB and mix.
- 4 Incubate at room temperature for 25 minutes.
- 5 Centrifuge.
- 6 Place on a magnetic stand until liquid is clear.
- 7 Remove and discard all supernatant.
- 8 Remove from the magnetic stand.
- 9 Add 200 µl EEW and mix.
- 10 Incubate at 50°C as follows.
  - ▶ [Plate] Place on the microheating system for 30 minutes.
  - ▶ [Tube] Place on the heat block for 30 minutes.
- 11 Place on a magnetic stand until liquid is clear.
- 12 Remove and discard all supernatant.
- 13 Remove from the magnetic stand.
- 14 Repeat steps 9–13 for a total of 2 washes.
- 15 Mix 28.5 µl EE1 and 1.5 µl HP3, and then vortex.
- 16 Add 23 µl elution premix and mix.
- 17 Incubate at room temperature for 2 minutes.
- 18 Centrifuge.
- 19 Place on a magnetic stand until liquid is clear.
- 20 Transfer 21 µl supernatant.
- 21 Add 4 µl ET2 and mix.
- 22 Centrifuge.

## Clean Up Captured Library

- 1 Add 45 µl SPB and mix.
- 2 Incubate at room temperature for 5 minutes.
- 3 Centrifuge.
- 4 Place on a magnetic stand until liquid is clear.
- 5 Remove and discard all supernatant.
- 6 Wash 2 times with 200 µl 80% EtOH.
- 7 Use a 20 µl pipette to remove residual EtOH.
- 8 Air-dry until dry.
- 9 Add 27.5 µl RSB and mix.
- 10 Remove from the magnetic stand.
- 11 Incubate at room temperature for 2 minutes.
- 12 Centrifuge.
- 13 Place on a magnetic stand until liquid is clear.
- 14 Transfer 25 µl supernatant.

## SAFE STOPPING POINT

If you are stopping, seal the plate or cap the tube and store at -25°C to -15°C for up to 7 days.

## Amplify Enriched Library

- 1 Add 5  $\mu$ l PPC.
- 2 Add 20  $\mu$ l EAM and mix.
- 3 Centrifuge.
- 4 Place on the thermal cycler and run the AMP10 program.

### SAFE STOPPING POINT

If you are stopping, seal the plate or cap the tube and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

## Clean Up Amplified Enriched Library

- 1 Centrifuge.
- 2 Transfer 50  $\mu$ l.
- 3 Add 50  $\mu$ l SPB and mix.
- 4 Incubate at room temperature for 5 minutes.
- 5 Centrifuge.
- 6 Place on a magnetic stand until liquid is clear.
- 7 Remove and discard all supernatant.
- 8 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 9 Use a 20  $\mu$ l pipette to remove residual EtOH.
- 10 Air-dry until dry.
- 11 Add 32  $\mu$ l RSB and mix.
- 12 Remove from the magnetic stand.
- 13 Incubate at room temperature for 2 minutes.
- 14 Centrifuge.
- 15 Place on a magnetic stand until liquid is clear.
- 16 Transfer 30  $\mu$ l supernatant.

### SAFE STOPPING POINT

If you are stopping, seal the plate or cap the tube and store at -25°C to -15°C for up to 7 days.

## Validate Enriched Libraries

- 1 Quantify using the Qubit dsDNA BR Assay Kit.
- 2 If the concentration is higher than the quantitative range for the High Sensitivity DNA chip, dilute the library 1:10 with RSB.
- 3 Run 1  $\mu$ l using a High Sensitivity DNA chip.

## Acronyms

Acronym	Definition
BLR	Blocker
CEX	Coding Exome Oligos
EAM	Enrichment Amplification Mix
EE1	Enrichment Elution Buffer 1
EEW	Enhanced Enrichment Wash Solution
EHB1	Enrichment Hybridization Buffer 1
EHB2	Enrichment Hybridization Buffer 2
ET2	Elute Target Buffer 2
HP3	2N NaOH
LAM	Library Amplification Mix
PPC	PCR Primer Cocktail
RSB	Resuspension Buffer
SMB	Streptavidin Magnetic Beads
SPB	Sample Purification Beads
ST2	Stop Tagment Buffer 2
TD	Tagment DNA Buffer
TDE2	Tagment DNA Enzyme 2