

DRAGEN TSO500 Analysis Software ICA Release Notes

V1.1.1

For TruSight Oncology 500 Assay

August 29, 2022



Introduction

These Release Notes detail the key features and known limitations to software components for the DRAGEN TSO500 v1.1.1 Analysis Software on Illumina Connected Analytics (ICA) platform.

This software is intended for use with the TruSight Oncology 500 Assay.

Software Version: 1.1.1

NEW FEATURES:

- Ability to run DRAGEN pipelines for TSO 500 analyses on ICA
- GUI enabled to run TSO500 analysis on ICA

KNOWN ISSUES:

- The sample sheet should not have blank rows between samples in the [Data] section, this may cause a run failure.
- Performance not verified using reads other than 2 x 101.
- The TSO500 RNA workflow is unstranded. Fusions or splice variants could involve antisense transcripts instead of the reported genes.
- Fusion caller may not always call fusions where breakpoint(s) are located in region(s) with high homology.
- The cloud workflow will fail if blank rows are present after the [Data] section.
- The Run ID field in the sample Information section for the SARJ files for some local DRAGEN analyses is displayed as "NA".
- The user interface in ICA allows the user to start an analysis without specifying the RUN folder path or FASTQ folder path. This action results in a failed analysis.
- An analysis run may be marked as 'Succeeded' although an analyzed sample may have failed.
 This may occur when a sample may have failed analysis, but the pipeline was executed
 successfully.

PRODUCT LIMITATIONS:

- The sample sheet must be configured as described in the User Guide.
- NovaSeq S4 analysis time for 192 samples may require more than 12 hours to complete.
- The values in the Run Metrics section will be listed as 'NA' if the analysis was started from FASTQs or if the analysis was started from BCLs but the InterOp files are missing or corrupted.
- Unmapped long insertions are not likely to occur on shorter indels because there is sufficient reference-matching sequence in the reads. Product claims only indels up to 25 base pairs.
- Complex variants are specifically output only for a specific region of the EGFR gene, component and phased variants would both be contained in the output.
- Incorrect calculation of variant allele frequency can occur in variants near the start and end of genomic reads, but there is a low probability of incorrect variant allele frequency in called variants due to sufficient variation in read start and end positions.



- Germline estimation uses latest publicly available population data and estimated to be representative of targeted population, the impact of rare germline mutations is expected to be limited.
- Germline estimation is difficult when tumor purity is >85% causing expected variant allele frequency for somatic and germline variants to converge.
- Poor quality wild type reads may align as chimeric and be miscalled during RNA analysis
- The Illumina Annotation Engine (aka Nirvana) may report incorrect protein (P-Dot) and transcript (C-Dot) changes in HGVS nomenclature for small variants located on a RefSeq transcript where an RNA-edit has occurred. Most known variants on these transcripts are unaffected. A list of affected Canonical RefSeq transcripts and Cosmic Variants from those transcripts can be found below. A full explanation of this product limitation can be found in PQN2020-1090. [1]

Affected Canonical RefSeq Transcripts

| Transcript ID | Gene Symbol |
|----------------|-------------|
| NM_002467.4 | MYC |
| NM_003224.5 | ARFRP1 |
| NM_004119.2 | FLT3 |
| NM_006904.6 | PRKDC |
| NM_198291.2 | SRC |
| NM_021960.4 | MCL1 |
| NM_001025366.2 | VEGFA |

Affected Cosmic Variants from Canonical RefSeq Transcripts



The list of affected variants is based on an analysis of COSMIC database version 92 variants located along the Canonical RefSeq Transcripts listed above [2]. New variants are regularly submitted to COSMIC, and this list of affected variants may change over time.

| Chr:Position | REF* | ALT** | Gene Symbol | Transcript ID | COSMIC_ID |
|----------------|------|--|----------------|--------------------|-------------------|
| chr1:150548890 | Α | АТСТА | MCL1 | NM_021960.4 | COSV57189 597 |
| chr6: 43738444 | С | Т | VEGFA | NM_0010253 66.2 | COSV10456 9261 |
| chr8:48805817 | G | GG | PRKDC | NM_006904.6 | COSV58041 377 |
| chr8:128748839 | GC | G | MYC | NM_002467.4 | COSV10438 8447 |
| chr8:128748840 | С | Α | MYC | NM_002467.4 | COSV10438 8806 |
| chr8:128748840 | С | G | MYC | NM_002467.4 | COSV10438 8204 |
| chr8:128748841 | Т | С | MYC | NM_002467.4 | COSV10438 8663 |
| chr13:28608094 | С | CACTITICCAAAAGCACCTGATCCTAGT ACCTTCCCAAACTCTAAATTTTCTCTTGG AAACTCCCATTTGAGATCATATTCATAT TCGTTCATC | FLT3 | NM_004119.2 | COSV54069 050 |
| chr13:28608124 | С | CTTCCCAAACTCTACTGTTGCGTTCATCA CTTTTCCAAAAGCACCTGATCCTAGTAC C | FLT3 | NM_004119.2 | COSV54044 227 |
| chr13:28608129 | С | CAAACTCAAAAGCACCTGATCCTAGTAC CTTCCC | FLT3 | NM_004119.2 | COSV54054 381 |
| chr13:28608129 | С | CAAACTCTAAATTTTCTCTTGGAAACTCC CATTATCCTAGTACCTTCCC | FLT3 | NM_004119.2 | COSV54043 729 |
| chr13:28608129 | С | CAAACTCTAAATTTTCTCTTGGAAACTCC CATTTTCCAAAAGCACCTGATCCTAGTA CCTTCCC | FLT3 | NM_004119.2 | COSV54075 746 |
| chr20:36030939 | G | GTGGCC | SRC | NM_198291.2 | COSV99050 886 |
| chr20:62331336 | С | сс | ARFRP1 | NM_003224.5 | COSV53926 174 |

^{*}Reference base(s)

[1] DRAGEN TSO500 uses the Canonical RefSeq transcript when annotating variants passed into the Combined Variant Output file. The Illumina Annotation Engine selects canonical transcripts based on the following rules:

- Order all overlapping transcripts by coding sequence length.
- Pick the longest transcript that has an associated Locus Reference Genome (LRG) sequence.

^{**}Alternate base(s)



- If no LRGs exist for the set of transcripts, pick the longest transcript that is coding.
- If there is a tie, pick the transcript with the smaller accession id number.

[2] Released 29 August 2022.

Release History

| Version | CN# | Author | Description of Change |
|---------|---------|-------------|-----------------------|
| 00 | 1073924 | Darryl Leon | Initial Release |