

# SurePlex Summary Protocol

## Reference Guide

Introduction .....	3
Workflow Diagram .....	5
Preparation and Materials .....	6
Preparation of Amplification Reaction Controls .....	7
Cell Lysis/Extraction .....	8
Pre-amplification .....	9
Amplification .....	10
Revision History .....	11
Technical Assistance .....	12



This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY, AND WILL VOID ANY WARRANTY APPLICABLE TO THE PRODUCT(S).

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2020 Illumina, Inc. All rights reserved.

All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see [www.illumina.com/company/legal.html](http://www.illumina.com/company/legal.html).

## Introduction

The SurePlex DNA Amplification System amplifies DNA from embryo biopsy samples (blastomere and trophectoderm biopsied cells) to quantities suitable for the VeriSeq PGS assay. This process creates highly representative samples that are ideal for pre-implantation genetic screening (PGS) during an *in vitro* fertilization (IVF) cycle.

SurePlex follows an easy-to-use single tube protocol. It is possible to complete the SurePlex protocol within a two and a half hour period.

## Recommendations

For quality results, adhere to the following recommendations.

### Sample Input Recommendations

All biopsy samples for VeriSeq PGS processing must be suspended in molecular grade 1x PBS (Phosphate Buffered Saline) and, if used, with a maximum concentration of 0.5% PVP (Polyvinylpyrrolidone) in a volume of 2.5  $\mu$ l. Validate before using any alternative sample collection volume, wash, and collection buffers.

Validate biopsy sample storage and shipment conditions before use. These factors may impact the quality and integrity of the samples. To maintain the integrity of the biopsy samples, store the biopsies at -65°C to -85°C and ship the biopsy tubes frozen. Store the biopsy tubes at -65°C to -85°C after collection. Complete the cell lysis to the PCR amplification steps within 14 days as described in the library preparation sections of this guide. Biopsy stability has not been evaluated under storage and shipping conditions other than these recommendations.

### Protocol Compliance

- ▶ Always use the current version of the SurePlex Summary Protocol Reference Guide, found on the Illumina Support website.
- ▶ Avoid contamination of samples and reagents during library preparation. Contamination of reagents will compromise your test results.
- ▶ Incubation times, incubation temperature, and pipetting volumes that differ from the specifications in this protocol will lead to suboptimal end products and compromise your test results.
- ▶ Adhere to all centrifuge and mixing steps in the protocol. Deviation from the recommended centrifugation and vortexing steps will lead to suboptimal end products and will compromise your test result.

### Warnings and Precautions

- ▶ Always use Personal Protective Equipment (PPE) during the laboratory steps.
- ▶ This kit is intended for use by qualified laboratory staff only. All laboratories are expected to follow Good Laboratory Practices (GLP) and have appropriate safety control measures in place.
- ▶ Check the documentation and safety data sheets (SDS) for guidelines on handling, preparing, and disposing kit components, intermediate mixtures, or wastes. For more information, see the SDS for this kit at [support.illumina.com/sds.html](https://support.illumina.com/sds.html).
- ▶ All components have an expiration date. Do not use kit components beyond the expiration date printed on component labels.

## Best Practices


- ▶ Conduct all pre-amplification activities (cell lysis, pre-amplification, and amplification set-up procedures) in a dedicated environment (Pre-PCR) physically separated from amplified genetic material.
- ▶ Use a sterile containment cabinet (vertical laminar flow cabinet) when handling embryo biopsies, during pre-amplification and amplification set up steps of the protocol. Follow applicable safety practices when handling embryo biopsies.
- ▶ Thoroughly clean the containment cabinet before starting the protocol.
- ▶ If recommended by the cabinet manufacturer, incubate tips, tubes, and pipettes with UV irradiation in the containment cabinet before starting the protocol.
- ▶ Make sure that you have a clean set of calibrated pipettes that are suitable for Pre-PCR activities, preferably exclusively dedicated to the SurePlex amplification protocol.

## Avoiding Contamination

- ▶ Pre-PCR environmental and sampling control measures (i.e., single-use of molecular grade reagents and consumables) are required to reduce the risk of sample contamination prior to input into the SurePlex reaction.
- ▶ Failure to follow Pre-PCR best practices during the reaction set-up may lead to contamination of the PCR products and may affect subsequent analysis and/or compromise your test results.

## Workflow Diagram

Figure 1 SurePlex Protocol Workflow Diagram



## Preparation and Materials

Table 1 Materials Required for Sample Preparation

Materials Required	Ordering Information
SurePlex DNA Amplification System	Illumina PR-40-415101-00
20x PBS (dilute to 1x with nuclease-free water)	Cell Signaling Technologies 9808
Genomic control DNA	Promega G1521
PCR tubes (0.2 ml, thin walled, flip cap) OR 96-well PCR plate and adhesive plate seals	Thermo Scientific AB-0620 Thermo Scientific AB-0600 (plate) AB-0558 (seals)
Microcentrifuge tube (1.5 ml, flip cap)	Sarstedt 72.690.001

Table 2 Starting Materials for Sample Preparation

Starting Materials
Single cell in 2.5 µl 1x PBS in 0.2 ml tubes or PCR plates
Positive control 1: 15 pg genomic DNA in 2.5 µl 1x PBS
Positive control 2: 60 pg genomic DNA in 2.5 µl 1x PBS
Negative control: 2.5 µl of 1x PBS (0.2 ml tube/PCR plate)
Negative control: 2.5 µl of collection buffer control (0.2 ml tube/PCR plate)

## SurePlex DNA Amplification System (50 reactions) - Kit Components

Table 3 SurePlex DNA Amplification System Kit Components

Component Name	Cap Color	Volume	Quantity	Prepare
Cell Extraction Buffer		300 µl	1	Thaw on ice
Extraction Enzyme Dilution Buffer		300 µl	1	Thaw on ice
Cell Extraction Enzyme		15 µl	1	Transfer to ice before use
SurePlex Pre-amp Buffer		300 µl	1	Thaw on ice
SurePlex Pre-amp Enzyme		15 µl	1	Transfer to ice before use
SurePlex Amplification Buffer		1.4 ml	1	Thaw on ice
SurePlex Amplification Enzyme		50 µl	1	Transfer to ice before use
Nuclease-free water	Clear	1 ml	2	Thaw on ice

## Preparation of Amplification Reaction Controls

### Estimated Time

Hands-on time: 30 minutes

### Consumables

Item	Quantity
1x PBS	805 $\mu$ l
Genomic DNA (100ng/ $\mu$ l)	5 $\mu$ l

### Steps

- 1 In a sterile containment cabinet (vertical laminar flow cabinet), label four sterile 1.5 ml microcentrifuge tubes and prepare their contents, as in [Table 4](#).

Table 4 Control Tubes for Sample Preparation

Control Tube	Description	Volume of 1x PBS	Volume of Female Genomic DNA /Stock Solution
1	1x PBS negative	100 $\mu$ l	None
2	2.5 ng/ $\mu$ l positive	195 $\mu$ l	5 $\mu$ l of 100 ng/ $\mu$ l stock
3	25 pg/ $\mu$ l positive	495 $\mu$ l	5 $\mu$ l of 2.5 ng/ $\mu$ l stock (from tube 2)
4	6.25 pg/ $\mu$ l positive	15 $\mu$ l	5 $\mu$ l of 25 pg/ $\mu$ l stock (from tube 3)

- 2 Add 1x PBS into each tube, adding the volumes shown in [Table 4](#).
- 3 Add 5  $\mu$ l of female genomic DNA to tube 2 and mix.
- 4 Dilute 5  $\mu$ l of tube 2 into tube 3 and mix.
- 5 Dilute 5  $\mu$ l of tube 3 into tube 4 and mix.
- 6 In the sterile containment cabinet, label three sterile 0.2 ml flat-capped PCR tubes as in [Table 5](#).

Table 5 PCR Tubes for Sample Preparation

PCR tube	Label on Tube Cap	Contents
1	15.6 pos	2.5 $\mu$ l of control tube 4 (6.25 pg/ $\mu$ l positive)
2	62.5 pos	2.5 $\mu$ l of control tube 3 (25 pg/ $\mu$ l positive)
3	Neg	2.5 $\mu$ l of control tube 1 (1x PBS)

- 7 Pipette 2.5  $\mu$ l of the control tubes into the PCR tubes, according to the contents listed in [Table 5](#).



#### NOTE

Use a clean tip for each tube and cap each tube immediately after transfer.

- 8 Store the control PCR tubes in a 96-well rack on ice until required.

## Cell Lysis/Extraction

### Estimated Time

Hands-on time: 30 minutes

Incubation time: 15 minutes

### Consumables

Item	Quantity
1x PBS	2.5 µl per sample
Cell Extraction Buffer (green cap)	2.5 µl per sample
Extraction Enzyme Dilution Buffer (purple cap)	4.8 µl per sample
Cell Extraction Enzyme (yellow cap)	0.2 µl per sample

### Steps

- 1 Collect samples in 2.5 µl of 1x PBS and 2.5 µl of collection buffer control. Line up sample tubes and control tubes in a 96-well rack on ice until required.
- 2 Centrifuge the tubes/plate at 200 × g for 3 minutes (at 4°C if possible).
- 3 Add 2.5 µl of Cell Extraction Buffer (green cap) to each sample (including controls) and store at 2°C to 8°C.
- 4 Following the specifications in [Table 6](#), prepare the Extraction Cocktail master mix on ice and mix well.

Table 6 Extraction Cocktail Master Mix

Extraction Cocktail	Cap Color	Volume per Single Sample	Volume per 5 Samples
Extraction Enzyme Dilution Buffer		4.8 µl	24 µl
Cell Extraction Enzyme		0.2 µl	1 µl
Total Volume		5 µl	25 µl

- 5 For every sample or control, add 5 µl of the freshly prepared Extraction Cocktail.
- 6 Briefly centrifuge samples to get all contents to the bottom of the tube.
- 7 Incubate samples in a PCR thermal cycler as in [Table 7](#).

Table 7 Thermal Cycler Program for Cell Lysis/Extraction

Number of Cycles	Temperature of Cycle	Incubation Time
1	75°C	10 min
1	95°C	4 min
1	22°C	Hold



#### CAUTION

Always use a heated lid on a PCR machine. Do not carry out PCR reactions under oil.



## Pre-amplification

### Estimated Time

Hands-on time: 10 minutes

Incubation time: 1 hour

### Consumables

Item	Quantity
SurePlex Pre-amp Buffer (red cap)	4.8 µl per sample
SurePlex Pre-amp Enzyme (white cap)	0.2 µl per sample

### Steps

- 1 Combine the Pre-amp Cocktail components as in [Table 8](#) and mix well.

Table 8 Pre-amp Cocktail

Pre-amp Cocktail	Cap Color	Volume per Single Sample	Volume per 5 Samples
SurePlex Pre-amp Buffer		4.8 µl	24 µl
SurePlex Pre-amp Enzyme		0.2 µl	1 µl
Total Volume		5 µl	25 µl

- 2 For each 10 µl sample prepared in *Cell Lysis/Extraction*, add 5 µl of SurePlex Pre-amp Cocktail. Briefly centrifuge.
- 3 Incubate samples according to the thermal cycler program in [Table 9](#).

Table 9 Thermal Cycler Program for Pre-amplification  
Reaction Products

Number of cycles	Temperature of cycle	Incubation Time
1	95°C	2 min
12	95°C	15 sec
	15°C	50 sec
	25°C	40 sec
	35°C	30 sec
	65°C	40 sec
	75°C	40 sec
1	4°C	Hold

- 4 Place the pre-amplification reaction products on ice.

## Amplification

### Estimated Time

Hands-on time: 10 minutes

Incubation time: 50 minutes



### Consumables

Item	Quantity
SurePlex Amplification Buffer (orange cap)	25 µl per sample
SurePlex Amplification Enzyme (blue cap)	0.8 µl per sample
Nuclease-free water (clear cap)	34.2 µl per sample

### Steps

- Combine the Amplification Cocktail components in [Table 10](#) and mix.

Table 10 Amplification Cocktail components

Amplification Cocktail	Cap Color	Volume per Single Sample	Volume per 5 Samples
SurePlex Amplification Buffer		25 µl	125 µl
SurePlex Amplification Enzyme		0.8 µl	4 µl
Nuclease-free water	Clear	34.2 µl	171 µl
<b>Total Volume</b>		<b>60 µl</b>	<b>300 µl</b>



#### NOTE

Sample amplification efficiency may be analyzed by using a real-time thermal cycler and adding SYBR Green I dye (Invitrogen, Catalog # S7563) at 0.125x final concentration in the Amplification Cocktail.

- Add 60 µl of the freshly prepared Amplification Cocktail to the 15 µl pre-amplification reaction product, cap the tube, and invert to mix. Centrifuge briefly.
- Amplify samples according to the thermal cycler program in [Table 11](#).

Table 11 Thermal Cycler Program for Amplification of Samples

Number of Cycles	Temperature of Cycle	Incubation Time
1	95°C	2 min
14	95°C	15 sec
	65°C	1 min
	75°C	1 min
1	4°C	Hold

- Keep amplified SurePlex products on ice before proceeding with the VeriSeq PGS assay protocol.

## SAFE STOPPING POINT

If you are stopping, leave the tubes on the thermal cycler overnight at 4°C or store at -25°C to -15°C for up to one month. The integrity of the PCR amplification products is affected by suboptimal storage conditions, the quality and grade of the plastic consumable used, and the nuclease-free control measured implemented..

## Revision History

Document	Date	Description of Change
Document # 15053626 v02	September 2020	<ul style="list-style-type: none"> <li>• Added Avoiding Contamination sub-section to Best Practices.</li> <li>• Updated workflow diagram.</li> <li>• Updated the lysis hold temperature program.</li> <li>• Updated the description of the safe stopping point after the amplification step.</li> <li>• Added the amplification hold temperature program.</li> </ul>
Document # 15053626 v01	February 2018	<ul style="list-style-type: none"> <li>• Converted document to new style and reference guide format.</li> <li>• Removed references to 24sure arrays.</li> <li>• Updated volumes in Kit Components table; added a Quantity column.</li> <li>• Revised Introduction and added Recommendations section.</li> <li>• Corrected minor errors and text edits.</li> </ul>
Document # 15053626 v00	July 2015	Initial release.

## Technical Assistance

For technical assistance, contact Illumina Technical Support.

Website: [www.illumina.com](http://www.illumina.com)  
 Email: [techsupport@illumina.com](mailto:techsupport@illumina.com)

### Illumina Technical Support Telephone Numbers

Region	Toll Free	International
Australia	+61 1800 775 688	
Austria	+43 800 006249	+43 1 9286540
Belgium	+32 800 77 160	+32 3 400 29 73
Canada	+1 800 809 4566	
China		+86 400 066 5835
Denmark	+45 80 82 01 83	+45 89 87 11 56
Finland	+358 800 918 363	+358 9 7479 0110
France	+33 8 05 10 21 93	+33 1 70 77 04 46
Germany	+49 800 101 4940	+49 89 3803 5677
Hong Kong, China	+852 800 960 230	
India	+91 8006500375	
Indonesia		0078036510048
Ireland	+353 1800 936608	+353 1 695 0506
Italy	+39 800 985513	+39 236003759
Japan	+81 0800 111 5011	
Malaysia	+60 1800 80 6789	
Netherlands	+31 800 022 2493	+31 20 713 2960
New Zealand	+64 800 451 650	
Norway	+47 800 16 836	+47 21 93 96 93
Philippines	+63 180016510798	
Singapore	1 800 5792 745	
South Korea	+82 80 234 5300	
Spain	+34 800 300 143	+34 911 899 417
Sweden	+46 2 00883979	+46 8 50619671
Switzerland	+41 800 200 442	+41 56 580 00 00
Taiwan, China	+886 8 06651752	
Thailand	+66 1800 011 304	
United Kingdom	+44 800 012 6019	+44 20 7305 7197
United States	+1 800 809 4566	+1 858 202 4566
Vietnam	+84 1206 5263	

**Safety data sheets (SDSs)**—Available on the Illumina website at [support.illumina.com/sds.html](http://support.illumina.com/sds.html).

**Product documentation**—Available for download from [support.illumina.com](http://support.illumina.com).



Illumina

5200 Illumina Way

San Diego, California 92122 U.S.A.

+1.800.809.ILMN (4566)

+1.858.202.4566 (outside North America)

[techsupport@illumina.com](mailto:techsupport@illumina.com)

[www.illumina.com](http://www.illumina.com)

For Research Use Only. Not for use in diagnostic procedures.

© 2020 Illumina, Inc. All rights reserved.

**illumina**<sup>®</sup>