

NextSeq 1000 and 2000

Denature and Dilute Libraries Guide

This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY, AND WILL VOID ANY WARRANTY APPLICABLE TO THE PRODUCT(S).

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2021 Illumina, Inc. All rights reserved.

All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.

Table of Contents

Overview	5
Loading Volume and Concentration	5
Recommended Loading Concentrations	5
Consumables and Equipment	6
Consumables	6
Equipment	6
Protocol	7
Prepare Reagents	7
Prepare a Fresh Dilution of NaOH	7
Prepare HT1	7
Dilute Libraries	7
Denature Libraries	9
Dilute Denatured Libraries	9
Dilute Denatured Libraries to Final Loading Concentration	9
Technical Assistance	10

Revision History

Document #	Date	Description of Change
1000000139235 v02	April 2021	Updated consumables.
1000000139235 v01	November 2020	Removed 10mM Tris HCL pH 8.5.
1000000139235 v00	October 2020	Initial release.

Overview

This guide explains how to manually denature and dilute prepared libraries for sequencing on the Illumina® NextSeq™ 1000 and NextSeq™ 2000 Sequencing Systems. Manual denature and dilution is recommended for low concentration libraries that are unable to meet the recommended loading concentration and volume stated in the *NextSeq 1000 and NextSeq 2000 Sequencing System Guide* (document # 1000000109376).

When initiating a run in the NextSeq 1000/2000 Control Software, make sure the Denature and Dilute On Board checkbox is not selected. See the *NextSeq 1000 and NextSeq 2000 Sequencing System Guide* (document # 1000000109376) for more information on run setup.

Loading Volume and Concentration

This procedure denatures and dilutes libraries to a final loading volume of 200 µl. Loading concentration can vary depending on library preparation and quantification methods.

Recommended Loading Concentrations

The optimal loading concentration depends on the library type and insert size. The following table provides the recommended final loading concentrations.

Library Type	Final Loading Concentration (pM)
AmpliSeq™ for Illumina Library PLUS	75
Illumina DNA Prep	75
Illumina DNA Prep with Enrichment	100
TruSeq DNA Nano 350	120
TruSeq DNA Nano 550	150
TruSeq Stranded mRNA	100
100 % PhiX	65
Illumina Stranded Total RNA with Ribo-Zero Plus	75

Library Type	Final Loading Concentration (pM)
Illumina Stranded mRNA Prep	75
Illumina DNA PCR-Free	100

Consumables and Equipment

Consumables

The following consumables are required to denature and dilute libraries.

Consumables	Supplier
1 N NaOH	General lab supplier
400 mM Tris-HCl, pH 8.0	General lab supplier
HT1 (Hybridization Buffer)	Illumina, catalog # 20015892 Illumina provided in the NextSeq 1000/2000 custom primers kits
Disposable gloves, powder-free	General lab supplier
Water, laboratory grade	General lab supplier
[Optional] PhiX Control v3	Illumina catalog # FC-110-3001

The following consumables for the denaturing and diluting libraries and PhiX are provided in the NextSeq 1000/2000 Reagents kits and Illumina library prep kits.

Consumables	Kit Name
NextSeq 1000/2000 RSB with Tween 20	Illumina provided in the NextSeq 1000/2000 P2 (v3) Reagents kit and the NextSeq 2000 P3 Reagents kit.
Microcentrifuge tube, 1.5 ml	VWR, catalog # 20170-038, or equivalent

Equipment

The following equipment is used to denature libraries that have been normalized.

Equipment	Supplier
Centrifuge	General lab supplier
Vortexer	General lab supplier

Protocol

Prepare Reagents

Prepare a Fresh Dilution of NaOH

Prepare a fresh dilution of 0.2 N NaOH to denature libraries for sequencing. Extra volume is prepared to prevent small pipetting errors from affecting the final NaOH concentration.

- Combine the following volumes in a microcentrifuge tube:
 - Laboratory-grade water (80 μ l)
 - Stock 1.0 N NaOH (20 μ l)

These volumes produce 100 μ l 0.2 N NaOH.

- Invert the tube several times to mix.

Prepare HT1

- Remove HT1 from -25°C to -15°C storage and thaw at room temperature.
- Store at 2°C to 8°C , or on ice, until you are ready to dilute denatured libraries.

Dilute Libraries

- Determine the required library concentration based on the desired final loading concentration.

Final Loading Concentration (pM)	Library Concentration (pM)
55	275
65	325
75	375
100	500
110	550

Final Loading Concentration (pM)	Library Concentration (pM)
120	600
130	650
140	700
150	750

2. Dilute libraries to the desired library concentration using RSB with Tween 20 to a final volume of 40 μ l.
3. [Optional] Add a PhiX Control.
 - a. Combine the following volumes in a microcentrifuge tube to prepare PhiX to match the library concentration.
 - 10 nM PhiX (2 μ l)
 - 10 mM RSB with Tween 20

Use the following table to determine the volume (μ l) of RSB with Tween 20 to use.

Initial PhiX Concentration (nM)	PhiX Volume (μ l)	Library Concentration (pM)	RSB With Tween 20 Volume (μ l)
10	2	275	70.7
10	2	325	59.5
10	2	375	51.3
10	2	500	38.0
10	2	550	34.4
10	2	600	31.3
10	2	650	28.8
10	2	700	26.6
10	2	750	24.7
10	2	1000	18.0

- b. Vortex briefly to mix, and then centrifuge at 280 x g for 1 minute.
- c. Add 1 μ l PhiX to 40 μ l nondenatured library diluted to the desired library concentration.

These volumes produce a ~2% PhiX spike-in. Actual percentage varies depending on library quality and quantity.

Denature Libraries

1. Add 10 μl 0.2 N NaOH to the tube of 40 μl nondenatured library and optional PhiX.
These volumes produce 50 μl denatured library, or 51 μl with PhiX.
2. Vortex briefly to mix, and then centrifuge at $280 \times g$ for 1 minute.
3. Incubate at room temperature for 8 minutes.

Dilute Denatured Libraries

1. Add 10 μl 400 mM Tris-HCl, pH 8.5 to the tube of denatured library.
These volumes produce 60 μl denatured library, or 61 μl with PhiX.
2. Vortex briefly to mix, and then centrifuge at $280 \times g$ for 1 minute.

Dilute Denatured Libraries to Final Loading Concentration

1. Add 140 μl prechilled HT1.
These volumes produce 200 μl final loading concentration, or 201 μl with PhiX.
2. Vortex to mix, and then centrifuge briefly.
3. Pierce the foil on the library reservoir.
4. Transfer the full volume of denatured and diluted library to the bottom of the cartridge reservoir provided with the NextSeq 1000/2000 P2 and P3 Reagents kits.

Technical Assistance

For technical assistance, contact Illumina Technical Support.

Website: www.illumina.com

Email: techsupport@illumina.com

Illumina Technical Support Telephone Numbers

Region	Toll Free	International
Australia	+61 1800 775 688	
Austria	+43 800 006249	+43 1 9286540
Belgium	+32 800 77 160	+32 3 400 29 73
Canada	+1 800 809 4566	
China		+86 400 066 5835
Denmark	+45 80 82 01 83	+45 89 87 11 56
Finland	+358 800 918 363	+358 9 7479 0110
France	+33 8 05 10 21 93	+33 1 70 77 04 46
Germany	+49 800 101 4940	+49 89 3803 5677
Hong Kong, China	+852 800 960 230	
India	+91 8006500375	
Indonesia		0078036510048
Ireland	+353 1800 936608	+353 1 695 0506
Italy	+39 800 985513	+39 236003759
Japan	+81 0800 111 5011	
Malaysia	+60 1800 80 6789	
Netherlands	+31 800 022 2493	+31 20 713 2960
New Zealand	+64 800 451 650	
Norway	+47 800 16 836	+47 21 93 96 93
Philippines	+63 180016510798	
Singapore	1 800 5792 745	
South Korea	+82 80 234 5300	

Region	Toll Free	International
Spain	+34 800 300 143	+34 911 899 417
Sweden	+46 2 00883979	+46 8 50619671
Switzerland	+41 800 200 442	+41 56 580 00 00
Taiwan, China	+886 8 06651752	
Thailand	+66 1800 011 304	
United Kingdom	+44 800 012 6019	+44 20 7305 7197
United States	+1 800 809 4566	+1 858 202 4566
Vietnam	+84 1206 5263	

Safety data sheets (SDSs)—Available on the Illumina website at support.illumina.com/sds.html.

Product documentation—Available for download from support.illumina.com.



Illumina

5200 Illumina Way

San Diego, California 92122 U.S.A.

+1.800.809.ILMN (4566)

+1.858.202.4566 (outside North America)

techsupport@illumina.com

www.illumina.com

For Research Use Only. Not for use in diagnostic procedures.

© 2021 Illumina, Inc. All rights reserved.

illumina[®]