HiSeq 1500 and HiSeq 2000 Systems Safety and Compliance Guide

For Research Use Only. Not for use in diagnostic procedures.

This guide provides important safety information pertaining to the installation, servicing, and operation of the HiSeq® 1500 and HiSeq® 2000 systems. This guide includes product compliance and regulatory statements. Read this document before performing any procedures on the system.

The country of origin and date of manufacture of the system are printed on the instrument label.

System Configuration

The HiSeq system is shipped with a dedicated control computer that is customized to system requirements. The computer is a dedicated subsystem of the instrument and is not intended for use or supported as a general-purpose computer.

For a list of components that are provided with the instrument, see the *HiSeq 2500, 1500, and 2000 Site Prep Guide* (document # 15006407).

Safety Considerations and Markings

This section identifies potential hazards associated with installing, servicing, and operating the instrument. Do not operate or interact with the instrument in a manner that exposes you to any of these dangers.

All of the hazards described herein can be avoided by following the standard operating procedures included in the *HiSeq 1500 System Guide* (document # 15035788) and the *HiSeq 2000 System Guide* (document # 15011190).

General Safety Warnings

Make sure that all personnel are trained in the correct operation of the instrument and any potential safety considerations.

Follow all operating instructions when working in areas marked with this label to minimize risk to personnel or the instrument.

Always exercise caution when the flow cell compartment door is open. The flow cell area contains components that could cause injury when the front door is open.

When the door is open, the X, Y, and Z motors are disabled by the instrument operating software.

Laser Safety Warning

The HiSeq 1500 and HiSeq 2000 systems are Class 1 laser products that each contain 2 Class 4 lasers.

Class 4 lasers present an eye hazard from direct and diffuse reflections. Avoid eye or skin exposure to direct or reflected Class 4 laser radiation. Class 4 lasers can cause combustion of flammable materials and produce serious skin burns and injury from direct exposure.

Do not operate the instrument with any of the panels removed. When the front door of the instrument is open, safety interlock switches block the laser beam. If you operate the instrument with any panels removed, you risk exposure to direct or reflected laser light.

Barcode Scanner Laser Precautions

A handheld barcode scanner that contains a Class 2 laser is included with the HiSeq 1500 and HiSeq 2000 systems. The barcode scanner is a Class 2 laser product. Do not stare into the visible-light beam of the barcode scanner.

Electrical Safety Warnings

Do not remove the outer panels from the instrument. There are no user-serviceable components inside. Operating the instrument with any of the panels removed creates potential exposure to line voltage and DC voltages.

The instrument is powered by 100–240 VAC operating at 50–60 Hz. Hazardous voltage sources are located behind the left back panel, but can be accessible if other panels are removed. Some voltage is present on the instrument even when the instrument is turned off. Operate the instrument with all panels intact to avoid electrical shock.

Power Specifications

Type	Specification
Line Voltage	100–240 VAC at 50–60 Hz
Power Consumption	Maximum 1500 Watts combined for the instrument, monitor, and work station

Electrical Connections

Connect the instrument to a grounded circuit capable of delivering at least:

- ▶ 10 Amps for a 100–110 Volt power source
- ▶ 6 Amps for a 220–240 Volt power source

For more information, see the *HiSeq 2500, 1500, and 2000 Site Prep Guide* (document # 15006407).

Protective Earth

The instrument has a connection to protective earth through the enclosure. The safety ground on the power cord returns protective earth to a safe reference. The protective earth connection on the power cord must be in good working condition when using this device.

Fuses

Only Illumina Field Service Engineers are qualified to replace the internal fuses. The power entry module includes 2 input fuses on the high-voltage input lines. These fuses are size 5x20 and are rated for 10 Amps, 250 VAC, Slo-Blo.

Hot Surface Safety Warning

Do not operate the instrument with any of the panels removed.

Do not touch the temperature station in the flow cell compartment. The heater used in this area is normally controlled between ambient room temperature (22°C) and 95°C. Exposure to temperatures at the upper end of this range can result in burns.

Heavy Object Safety Warning

The instrument is heavy and can cause serious injury if dropped or mishandled.

Only personnel authorized by Illumina can uncrate, install, or move the instrument. Each time the instrument is moved, install the shipping brackets to prevent stage assembly damage. Significant risks to optical and mechanical alignment can occur. The instrument must be properly repositioned and recalibrated.

Uncrating, Installing, and Moving the Instrument

Only personnel authorized by Illumina can uncrate, install, or move the instrument. If you have to relocate the instrument, contact your Illumina representative.

Environmental Considerations

Element	Specification
Temperature	Maintain a lab temperature of 19°C to 25°C (22°C ±3°C). This temperature is the operating temperature of the instrument. During a run, do not allow the ambient temperature to vary more than ±2°C.
Humidity	Maintain a noncondensing relative humidity between 20–80%.
Elevation	Locate the instrument at an altitude below 2000 meters (6500 feet).
Air Quality	Operate the instrument in a Pollution Degree II environment or better. A Pollution Degree II environment is defined as an environment that normally includes only nonconductive pollutants.
Ventilation	Consult your facilities department for ventilation requirements based on the instrument heat output specifications.

Product Certifications and Compliance

The HiSeq 1500 and HiSeq 2000 are certified to the following standards:

- ▶ UL STD 61010-1
- ▶ CSA STD C22.2 No. 61010-1
- ▶ IEC/EN 61010-1
- ▶ IEC/EN 61326-1

The HiSeq 1500 and HiSeq 2000 comply with the following directives:

- ▶ Low Voltage Directive 2006/95/EC
- ▶ EMC Voltage Directive 2004/108/EC

Compliance and Regulatory Markings

The instrument is labeled with the following compliance and regulatory markings.

This label assures that the product is compliant with safety standards and has been tested and certified by a third-party organization.

This label assures that the product meets the essential requirements of all relevant EU directives.

This label assures that the product complies with the Environmental Protection Use Period - 10 years.

Restriction of Hazardous Substances (RoHS)

This label indicates that the instrument meets the WEEE Directive for waste.

Return the instrument to Illumina for disposal.

Product Compliance and Regulatory Statements

FCC Compliance

This device complies with Part 15 of the FCC Rules. Operation is subject to the following 2 conditions:

- 1 This device may not cause harmful interference.
- 2 This device must accept any interference received, including interference that may cause undesired operation.

CAUTION

Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instrumentation manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case users will be required to correct the interference at their own expense.

Shielded Cables

Shielded cables must be used with this unit to ensure compliance with the Class A FCC limits.

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Notes

Revision History

Document	Date	Description of Change
Document # 1000000000665 v00	February 2016	Initial release.

Copyright and Trademarks

© 2016 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiniSeq, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.