NeoPrep Library Prep System Safety and Compliance Guide

For Research Use Only. Not for use in diagnostic procedures.

This guide provides important safety information pertaining to the installation, servicing, and operation of the Illumina® NeoPrep™ Library Prep System. This guide includes product compliance and regulatory statements. Read this document before performing any procedures on the system.

The country of origin and date of manufacture of the system are printed on the instrument label.

Safety Considerations and Markings

This section identifies potential hazards associated with installing, servicing, and operating the instrument. Do not operate or interact with the instrument in a manner that exposes you to any of these dangers.

All of the hazards described herein can be avoided by following the standard operating procedures included in the *NeoPrep Library Prep System Guide* (document # 15049720).

General Safety Warnings

Make sure that all personnel are trained in the correct operation of the instrument and any potential safety considerations.

Follow all operating instructions when working in areas marked with this label to minimize risk to personnel or the instrument.

Electrical Safety Warnings

Do not remove the outer panels from the instrument. There are no user-serviceable components inside. Operating the instrument with any of the panels removed creates potential exposure to line voltage and DC voltages.

The instrument is powered by 100–240 VAC operating at 50–60 Hz. Hazardous voltage sources are located on the rear bottom left, but can be accessible if other panels are removed. Some voltage is present on the instrument even when the instrument is turned off. Operate the instrument with all panels intact to avoid electrical shock.

Power Specifications

- ▶ Line Voltage: 100–240 Volts AC @ 50/60 Hz
- ▶ Power Consumption: 220 Watts

Electrical Connections

Connect the instrument to a grounded circuit capable of delivering at least:

- ▶ 15 Amps for a 100–120 Volt power source
- ▶ 10 Amps for a 220–240 Volt power source

For more information, see the *NeoPrep System Site Prep Guide* (document # 15050812).

Protective Earth

The instrument has a connection to protective earth through the enclosure. The safety ground on the power cord returns protective earth to a safe reference. The protective earth connection on the power cord must be in good working condition when using this device.

Fuses

Only Illumina Field Service Engineers are qualified to replace the internal fuses. The power entry module includes 2 input fuses on the high-voltage input lines. These fuses are rated for 7 Amps, 250 VAC, Slo-Blo.

Hot Surface Safety Warning

Do not operate the instrument with any of the panels removed.

When the instrument is performing a run, do not touch the surfaces of the integrated library prep module that are accessible through the rear of the library card compartment. Exposure to hot surfaces could result in burns.

Heavy Object Safety Warning

The instrument weighs approximately 21 kg (46 lb) and can cause serious injury if dropped or mishandled. A 2 person lift is required.

Uncrating, Installing, and Moving the Instrument

Only personnel authorized by Illumina can uncrate, install, or move the instrument. If you have to relocate the instrument, contact your Illumina representative.

Environmental Considerations

Element	Specification
Temperature	Maintain a lab temperature of 19°C to 25°C (22°C ±3°C). This temperature is the operating temperature of the instrument.
Humidity	Maintain a noncondensing relative humidity between 20–80%.
Elevation	Locate the instrument at an altitude below 2000 meters (6500 feet).
Air Quality	Operate the instrument in a Pollution Degree II environment or better. A Pollution Degree II environment is defined as an environment that normally includes only nonconductive pollutants.
Ventilation	Consult your facilities department for ventilation requirements based on the instrument heat output specifications.

Product Certifications and Compliance

The NeoPrep System is certified to the following standards:

- ▶ CE marked
- ▶ UL STD 61010-1
- ▶ CSA STD C22.2 No 61010-1
- ▶ EN/IEC 61010-1
- ▶ EN/IEC 61326-1
- ▶ EN/IEC 61326-2-6

The NeoPrep System complies with the following directives:

- ▶ Low Voltage Directive 2006/95/EEC
- ▶ EMC Directive 2004/108/EC

Compliance and Regulatory Markings

The NeoPrep System is labeled with the following compliance and regulatory markings.

This label assures that the product is tested and certified by TUV Rheinland, a Nationally Recognized Testing Laboratory (NRTL).

This label assures that the product meets the essential requirements of all relevant EU directives.

This label assures that the product complies with the Environmental Protection Use Period - 10 years.

Restriction of Hazardous Substances (RoHS)

This label indicates that the instrument meets the WEEE Directive for waste.

Return the instrument to Illumina for disposal.

Product Compliance and Regulatory Statements

FCC Compliance

This device complies with Part 15 of the FCC Rules. Operation is subject to the following 2 conditions:

- 1 This device may not cause harmful interference.
- 2 This device must accept any interference received, including interference that may cause undesired operation.

CAUTION

Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instrumentation manual, may cause harmful interference to radio communications.

Operation of this equipment in a residential area is likely to cause harmful interference in which case users will be required to correct the interference at their own expense.

Shielded Cables

Shielded cables must be used with this unit to ensure compliance with the Class A FCC limits.

IC Compliance

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

This device complies with Industry Canada licence-exempt RSS standards. Operation is subject to the following two conditions:

- 1 This device may not cause interference.
- 2 This device must accept any interference, including interference that may cause undesired operation of the device.

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Revision History

Document	Date	Description of Change
Material # 20000942 Document # 15050811 v01	October 2015	Added a revision history Added translations for French, German, Italian, Korean, Russian, and Spanish. Added the Korea compliance statement.
Part # 15050811 Rev. B	March 2015	Corrected humidity operating conditions to 20–80% relative humidity (noncondensing).
Part # 15050811 Rev. A	March 2015	Initial release.

Copyright and Trademarks

© 2015 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.