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
A.

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B.

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Illumina
Complete Long
Read Prep with
Enrichment

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[Illumina Complete Long Read Prep with Enrichment, Human](#)

[Human Comprehensive Panel](#)

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表 3: Human Comprehensive Panel 性能

Human Comprehensive Panel 目标区域的准确率 ^a			
	Illumina Complete Long Read Prep with Enrichment 和 Human Comprehensive Panel ^b	标准短读长 WGS ^c	
F1 评分 (SNV)	99.54%	99.04%	
F1 评分 (插入缺失)	99.47%	99.15%	
F1 评分 (SV)	80.43%	60.18%	
全基因组、所有基准区域的准确率 ^d			
F1 评分 (SNV+ 插入缺失)	Illumina Complete Long Read Prep with Enrichment 和 Human Comprehensive Panel ^b	Illumina Complete Long Read Prep, Human ^e	标准短读长 WGS ^c
	99.87%	99.90%	99.84%

a. 仅常染色体。SNV, 单核苷酸位点变异; SV, 结构变异。
b. Illumina Complete Long Read Prep with Enrichment, Human Comprehensive Panel + Illumina DNA PCR-Free Prep + NovaSeq 6000 基因测序仪 + DRAGEN v4.2。
c. Illumina DNA PCR-Free Prep + NovaSeq 6000 基因测序仪 + DRAGEN v4.2。
d. 使用 PrecisionFDA Truth Challenge V2 基准数据² 和 DRAGEN v4 的内部存档数据, 与所有参赛解决方案进行比较。
e. Illumina Complete Long Read Prep, Human + Illumina DNA PCR-Free Prep + NovaSeq 6000 基因测序仪 + DRAGEN v4.2。

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1. Bekritsky Ma, Bekritsky Ma, Colombo C, Eberle Ma. Identifying genomic regions with high quality single nucleotide variant calling. illumina.com/science/genomics-research/articles/identifying-genomic-regions-with-high-quality-singlenucleotide-.html. Published 2021. Accessed August 30, 2023.
2. PrecisionFDA. Truth Challenge V2: Calling Variants from Short and Long Reads in Difficult-to-Map Regions. [Precision.fda.gov/challenges/10](https://precision.fda.gov/challenges/10). Accessed October 2, 2023.
3. Roessler K. Illumina Complete Long Reads software analysis workflow for human WGS. illumina.com/science/genomicsresearch/articles/complete-long-read-software-analysis.html. Published 2023. Accessed September 22, 2023.

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产品	货号
Illumina Complete Long Read Prep with Enrichment, Human (24 样本)	20113832
Illumina Complete Long Read Prep with Enrichment, Human (96 样本)	20113833
Illumina Complete Long Read Prep with Enrichment, Human Comprehensive Panel (24 样本)	20113834
Illumina Complete Long Read Prep with Enrichment, Human Comprehensive Panel (96 样本)	20113835
Illumina Human Comprehensive Panel (24 样本)	20113836
Illumina Human Comprehensive Panel (96 样本)	20113837
Illumina Unique Dual Indexes, LT (48 标签, 48 样本)	20098166
Illumina DNA/RNA UD Indexes Set A, Tagmentation (96 标签, 96 样本)	20091654
Illumina DNA/RNA UD Indexes Set B, Tagmentation (96 标签, 96 样本)	20091656
Illumina DNA/RNA UD Indexes Set C, Tagmentation (96 标签, 96 样本)	20091658
Illumina DNA/RNA UD Indexes Set D, Tagmentation (96 标签, 96 样本)	20091660
Illumina Analytics - 1 iCredit	20042038

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© 2024 Illumina, Inc. 保留所有权利。所有商标均为因美纳公司或其各自所有者的财产。