

# Methods for RNA sequencing

Illumina solutions for profiling  
RNA, from targeted panels to  
the whole transcriptome



Table of contents

Introduction

Key benefits ..... 4

Overview of the RNA-Seq workflow ..... 4

Method 1: mRNA-Seq

Relevant applications ..... 8

Step-by-step overview ..... 9

Method 2: Targeted RNA-Seq

Relevant applications ..... 12

Step-by-step overview ..... 13

Method 3: Total RNA-Seq

Relevant applications ..... 16

Step-by-step overview ..... 17

Support that never stops

Trusted technology, trusted partner ..... 20

Learn more..... 20

Ordering information

Library preparation..... 21

Sequencing systems..... 22

Data analysis..... 22

Appendix

References..... 23


Introduction

As researchers gain deeper insights into the genome, they learn how factors beyond the genetic sequence influence development, disease states, and more. To truly understand what is occurring, scientists are turning towards a multiomics approach, looking at the transcriptome, epigenome, and proteome, in addition to the genome, to find comprehensive answers. In this guide, we will focus on next-generation sequencing (NGS) workflows for studying the transcriptome.

NGS-based RNA sequencing (RNA-Seq) offers a highly sensitive method for analyzing gene expression across the transcriptome, illuminating changes and characterizing various RNA forms with no prior knowledge of the transcriptome required. This technique uncovers transcript isoforms, gene fusions, and single nucleotide variants (SNVs) in a single experiment.<sup>1-3</sup> Bulk RNA-Seq is a well-established method that measures average RNA expression in cell populations. Ideal for newcomers to NGS, bulk RNA-Seq is increasingly used in translational and clinical cancer research, benefiting from its simplicity and the depth of information RNA-Seq provides.<sup>4</sup>

Overview of the RNA-Seq workflow

Illumina RNA-Seq workflows integrate RNA extraction, library preparation, sequencing, and data analysis to support transcriptome studies (Figure 1).



Key benefits:

**Comprehensive analysis:** RNA-Seq offers several advantages over other RNA analysis methods, such as qPCR and gene expression arrays, including:

- Covering a wider dynamic range than gene expression arrays, resulting in greater sensitivity and accuracy<sup>5</sup>
- Capturing both known and novel features, enabling analysis of the transcriptome without a reference<sup>2</sup>
- Offering a rich view of the transcriptome beyond expression profiling, detecting alternative splice sites, gene fusions, and allele-specific expression

STEP 1  
Nucleic acid extraction

Many commercial DNA and RNA extraction and purification kits are available. Choose a kit or method that produces the highest quality nucleic acids possible for your specific targets and from your specific sample types. Illumina has tested a variety of commercially available kits, but other kits or methods work better for your targets and samples.

STEP 2  
Library preparation

Following nucleic acid isolation, scientists can prepare nucleic acid libraries—a collection of similarly sized fragments that have known oligonucleotide adapter sequences attached to the 5' and 3' ends of the strands that are loaded onto an instrument for sequencing. The Illumina RNA library prep portfolio offers a range of solutions to support various sequencing methods and sample types, including formalin-fixed paraffin-embedded (FFPE) tissue (Table 1). Illumina RNA library prep solutions offer flexibility, scalability, and performance with a rapid, automation-friendly workflow option to prepare sequencing-ready libraries in a single day. The RNA library preparation products mentioned in this guide are representative of the portfolio and are not meant to be a comprehensive list of all available options. View the full RNA library prep portfolio.

Table 1: A comparison of select Illumina RNA preparation kits<sup>a</sup>

Parameter	Illumina Stranded mRNA Prep	Illumina RNA Prep with Enrichment	Illumina Stranded Total RNA Prep with Ribo-Zero Plus
Area of interest	Coding transcriptome with high-quality samples	Coding transcriptome with FFPE samples	Noncoding and coding RNA
Description	Quantify gene expression, identify known and novel isoforms, detect gene fusions, and measure allele-specific expression	Analyze expression in a focused set of genes of interest, obtain quantitative expression information, detect small variants and gene fusions	Measure genes and transcripts while detecting known and novel features; targeted hybridization removes abundant rRNA to focus on high value portions of the transcriptome
Mechanism of action	PolyA selection upfront of ligation-based RNA library prep	Tagmentation-based RNA library prep followed by hybrid-capture enrichment	rRNA depletion ahead of ligation-based RNA library prep
Input	25–1000 ng high-quality RNA	10 ng high-quality RNA, 20 ng FFPE RNA	1 ng high-quality RNA, 10 ng FFPE RNA
Hands-on time	< 3 hr	< 2 hr	< 3 hr
Total assay time	6.5 hr	< 9 hr	~ 7 hr
Automation capability	Yes	Yes	Yes
a. More RNA library preparation kits are available at <a href="https://illumina.com/techniques/sequencing/ngs-library-prep/rna.html">illumina.com/techniques/sequencing/ngs-library-prep/rna.html</a> .			



STEP 3  
Sequencing

After libraries are prepared, they are ready for sequencing. Regardless of your research question, flexible Illumina sequencing systems can help you find answers using simple push-button workflows ([Table 2](#)). For studies focused on a small amount of information, ie, prokaryotic species or small targeted RNA oncology panels, researchers can use a benchtop sequencing system, such as the MiSeq™ i100 Series or NextSeq™ 1000 and NextSeq 2000 Systems. For large-scale studies, researchers can use high-throughput instruments like the NovaSeq™ 6000 System and NovaSeq X Series, and multiplex up to 384 samples with unique dual indexes.

Table 2: Select examples of sequencing systems



System	MiSeq i100 Series	NextSeq 1000 and NextSeq 2000 Systems	NovaSeq 6000 System	NovaSeq X Series
System overview	Simplified operations and fast, flexible sequencing	Expansive application breadth and proven performance	Immense discovery power for deeper insights	Extraordinary throughput and transformative economics
Output range	1.5–30 Gb	10–540 Gb <sup>a</sup>	65 Gb–6 Tb	165 Gb–16 Tb
Single reads per flow cell	5–100M	100M–1.8B	800M–10B	1.6–26B
Maximum read length	2 × 500 bp	2 × 300 bp	2 × 250 bp	2 × 150 bp
a. Maximum specifications based on a P4 flow cell run; P4 flow cells are available for the NextSeq 2000 System only.				

STEP 4  
Data analysis

RNA-Seq data can be easily and securely transferred, stored, and analyzed in the Illumina cloud-computing platforms: BaseSpace™ Sequence Hub and Illumina Connected Analytics. For users without bioinformatics expertise, BaseSpace Sequence Hub is recommended for its intuitive interface, easy run setup and monitoring, and simplified push-button secondary analysis. For more advanced users, Illumina Connected Analytics supports customization with highly configurable, scalable analysis. Both platforms offer in-cloud access to DRAGEN™ secondary analysis pipelines for accurate and efficient analysis of RNA-Seq data.

STEP 5  
Insights

Outputs from secondary analysis pipelines can be ingested into reporting and exploration software, such as:

- [Correlation Engine](#)—analyze private omics data with highly curated public data to help put data into biological context
- [Illumina Connected Insights](#)—streamline variant interpretation and reporting for somatic oncology research applications
- [Partek™ Flow™ software](#)—perform statistical analysis and explore data with interactive, customizable, and publication-ready visualizations



## Method 1: mRNA-Seq

Messenger RNA sequencing (mRNA-Seq) sensitively and accurately quantifies gene expression, identifies known and novel isoforms in the coding transcriptome, and measures allele-specific expression. Protein-coding genes with polyA tails are selected ahead of library preparation in this method. Using mRNA-Seq to study the coding transcriptome, researchers can focus on a smaller, more manageable portion of the transcriptome that will provide information relevant to their area of interest.

### Relevant applications

#### Gene expression profiling for disease research

To understand normal cell development and disease mechanisms, researchers frequently investigate differential expression during development, in specific tissues, or in response to varying conditions. mRNA-Seq shows exceptional performance in profiling genes with low expression levels. It is being used to assess gene expression profiles for the study of complex diseases and laying the groundwork for advances in precision medicine by identifying potentially therapeutic biomarkers.<sup>6</sup>

Data from mRNA-Seq experiments can offer insight into the gene networks and pathways involved in complex disease and cell biology mechanisms.<sup>7,8</sup> For example, transcriptomic analysis is helping researchers compare brain regions with different pathology to identify meaningful gene expression changes in Alzheimer's disease (AD).<sup>7,9</sup> Differential expression profiling is also revealing the pathogenesis of heart failure and identifying gene signatures to detect heart disease.<sup>10-14</sup>

#### Biomarker profiling for drug development


mRNA-Seq is increasingly being utilized to discover and profile RNA-based drug response biomarkers with the aim of improving the efficiency and success rate of the drug development process. While several technologies have been used for this application, the capabilities of mRNA-Seq promise to be of particular benefit.<sup>15-17</sup>

Differential expression profiles and gene fusions detectible through RNA analysis have been shown to associate with a range of response characteristics, including efficacy, the incidence of adverse effects, pharmacodynamics, and other attributes.<sup>18-21</sup> Such biomarkers have therefore become an invaluable tool in multiple components of the development process, such as informing the interpretation of clinical trial data, allowing more efficient stratification of trial cohorts, and identifying neoantigen candidates for immunotherapy.<sup>22,23</sup> RNA-based biomarkers may also provide a foundation for the development of companion diagnostics, including for compounds that have previously failed clinical trials.



### Step-by-step overview

mRNA-Seq has five basic steps: RNA extraction, library preparation, sequencing, analysis, and insights (Figure 2).



#### STEP 1

##### RNA extraction

For mRNA-Seq, it is crucial to use high-quality mRNA as input. There are several commercial options for mRNA isolation, depending on the sample type; the QIAGEN RNeasy for RNA kit is recommended for use with Illumina mRNA-Seq workflows. Effectively evaluating mRNA quality is a critical step in successful RNA sequencing and can be achieved by measuring mean mRNA fragment size.

#### STEP 2

##### Library preparation

Illumina Stranded mRNA Prep delivers accurate, unbiased detection of the coding transcriptome with precise measurement of strand information (Table 3).

Table 3: Illumina targeted amplicon sequencing analysis software

Feature	Specification
RNA input amount	25–1000 ng high-quality RNA
Total assay time	6.5 hr
Hands-on time	< 3 hr



STEP 3  
Sequencing

Several sequencing systems can be used for mRNA-Seq. The one chosen depends on several factors, including the application, study size, throughput requirements, and more (Table 4).

Table 4: Experimental parameters for performing mRNA-Seq on different sequencing systems

System	Flow cell	Single reads per flow cell	No. of samples per flow cell <sup>a</sup>	Recommended read length
MiSeq i100 Series	50M	50M	2	2 × 75 bp
	100M	100M	4	
NextSeq 1000 and NextSeq 2000 Systems	P1	100M	4	2 × 75 bp
	P2	400M	16	
	P3 <sup>b</sup>	1.2B	48	
	P4 <sup>b</sup>	1.8B	72	
NovaSeq 6000 System	SP	1.6B	64	2 × 75 bp
	S1	3.2B	128	
	S2	8.2B	328	
	S4	8–10B	384	
NovaSeq X Series	10B	10B	384	2 × 75 bp
	25B	26B	384	

a. Based on 25M reads per sample. Sufficient gene expression and most use cases of fusion calling.

b. Flow cells are only available on the NextSeq 2000 System.

STEP 4  
Data analysis

For secondary analysis of mRNA-Seq data, Illumina recommends using the DRAGEN RNA pipeline, which performs RNA quantification, gene fusion detection, and small variant calling in one integrated workflow (Table 5). This pipeline is available on-premises with a DRAGEN Server, on the cloud-based BaseSpace Sequence Hub and Illumina Connected Analytics platforms, and onboard select sequencing systems, including the NextSeq 1000 and NextSeq 2000 Systems and NovaSeq X Series.

STEP 5  
Insights

After secondary data analysis, results can be transferred to Correlation Engine to understand the biological effects of gene expression changes (Table 5). Correlation Engine contains knowledge-based gene sets and results from thousands of public studies that inform biological interpretation. Connect differential gene expression data from RNA-Seq experiments with disease associations or visualize correlated genes. Gain further insights with Partek Flow software using robust statistical algorithms and information-rich visualizations. Additionally, Illumina Connected Insights can be used to streamline interpretation and reporting for research applications.

Table 5: Illumina mRNA-Seq analysis software

Pipeline	Application	Input	Access point
<a href="#">DRAGEN RNA</a>	Offers multiple operating modes, including reference-only alignment and annotation-assisted alignment with gene fusion detection. The gene fusion module uses the DRAGEN RNA spliced aligner to perform split-read analysis on supplementary (chimeric) alignments to detect potential breakpoints.	FASTQ  Optional: Custom reference Gene Annotation File (GTF, GFF, or GFF3)	<ul style="list-style-type: none"><li>DRAGEN Server</li><li>BaseSpace Sequence Hub</li><li>Illumina Connected Analytics</li><li>Onboard MiSeq i100 Series, NextSeq 1000/2000 Systems, NovaSeq X Series</li></ul>
<a href="#">RNA-Seq interpretation</a>	Processes sequencing data from mRNA to estimate transcript abundance and identify differentially expressed transcripts across samples	csv, txt, xlsx files	<ul style="list-style-type: none"><li>Correlation Engine</li></ul>
<a href="#">RNA-Seq visualization and statistical analysis</a>	Provides differential analysis, clustering, and data exploration plots	tsv, csv, txt, gz, FASTQ, BAM	<ul style="list-style-type: none"><li>Partek Flow software</li></ul>
<a href="#">Illumina Connected Insights</a>	Supports streamlined interpretation and reporting from DRAGEN software for oncology research applications	VCF	<ul style="list-style-type: none"><li>Illumina Connected Insights with automated ingestion of VCF files</li></ul>

Method 2: Targeted RNA-Seq

Targeted RNA-Seq is a highly accurate method for selecting and sequencing specific transcripts of interest that offers both quantitative and qualitative information. Targeted RNA-Seq can be achieved via either enrichment or amplicon-based approaches, both of which enable gene expression analysis in a focused set of genes of interest. Here, we focus on enrichment assays, which provide the ability to detect both known and novel gene fusion partners in many sample types, including FFPE tissue.

Relevant applications

Variant detection for cancer research

Targeted RNA-Seq is a critical tool for direct measurement of the functional consequences of mutations. Despite the average cancer containing about 46 mutations, only five to eight are necessary for initiation.<sup>24</sup> Genomic profiling alone is insufficient to differentiate these driver mutations from passenger mutations, or those that do not influence cancer initiation or progression. Measurement of gene expression patterns and mutation consequences using targeted RNA-Seq enables large-scale, unbiased differentiation of factors crucial for cancer progression, resulting in more thorough and accurate cancer modeling.


Detection and characterization of respiratory pathogens

Targeted RNA-Seq provides an effective method for rapid and accurate identification of respiratory pathogens. Combining Illumina RNA Prep with Enrichment with target-specific probe panels enables targeted sequencing of different subsets of respiratory pathogens. For example, the Respiratory Virus Oligo Panel v2 targets SARS-CoV-2 and other common respiratory viruses in a single assay.<sup>25</sup> Alternatively, Illumina offers the Respiratory Pathogen ID/AMR Panel, which targets ~280 respiratory pathogens, including viruses, bacteria, and fungi, and associated antimicrobial resistance (AMR) markers.<sup>26</sup>



Step-by-step overview

Targeted RNA-Seq has five basic steps: RNA extraction, library preparation and enrichment, sequencing, data analysis, and insights (Figure 3).



STEP 1  
RNA extraction

For targeted RNA-Seq, high-quality intact RNA or RNA isolated from FFPE samples can be used as input. There are several commercial options for RNA isolation, depending on the sample type; the QIAGEN RNeasy for RNA kit is recommended for use with intact samples and QIAGEN RNeasy FFPE for use with FFPE tissue. Effectively evaluating RNA quality is a critical step in successful RNA-Seq and can be achieved by measuring mean RNA fragment size before library prep.

STEP 2  
Library preparation

Illumina RNA Prep with Enrichment provides a fast, integrated workflow for producing enriched and indexed sequencing libraries (Table 6). The kit enables rapid, targeted interrogation of an expansive number of target genes with exceptional capture efficiency and coverage uniformity using various prebuilt probes or custom probes for maximum flexibility.<sup>27</sup>

Table 6: Illumina RNA Prep with Enrichment specifications

Feature	Specification
RNA input amount	10 ng high-quality RNA, 20 ng FFPE RNA
Total assay time	< 9 hr
Hands-on time	< 2 hr

STEP 3  
Sequencing

There are many panels that can be used to select the targets of interest. In this example workflow, we focus on using the Illumina Exome Panel to analyze the coding transcriptome. The sequencing system used depends on several factors, including the application, study size, throughput requirements, and more ([Table 7](#)).

Table 7: Experimental parameters for performing targeted RNA-Seq on different sequencing systems

System	Flow cell	Single reads per flow cell	No. of samples per flow cell <sup>a</sup>	Recommended read length
MiSeq i100 Plus System	25M	25M	8	2 × 100 bp
	50M	50M	16	
	100M	100M	33	
NextSeq 1000 and NextSeq 2000 Systems	P1	100M	33	2 × 100 bp
	P2	400M	133	
a. Based on 3M clusters per sample using the TruSight RNA Pan-Cancer Panel.				

STEP 4  
Data analysis

Illumina recommends using the DRAGEN RNA pipeline or the Enrichment panel-specific analysis workflow, ie, DRAGEN Microbial Enrichment Plus ([Table 8](#)). These pipelines are available on-premises with a DRAGEN Server, on cloud-based platforms, including BaseSpace Sequence Hub and Illumina Connected Analytics, and onboard select sequencing systems.

STEP 5  
Insights

After secondary analysis, further explore data with Partek Flow software using statistical analysis and interactive visualizations. For somatic oncology applications using targeted RNA-Seq, Illumina Connected Insights enables QC, annotation, interpretation, curation of SNVs, splice variants, and fusion variants from RNA (expression analysis coming in future) with subsequent report generation. RNA may be assessed together with DNA or sequentially, and with either option, results may be merged into a single research report if desired.

Table 8: Illumina mRNA-Seq analysis software

Pipeline	Application	Input	Access point
<a href="#">DRAGEN RNA</a>	Offers multiple operating modes, including reference-only alignment and annotation-assisted alignment with gene fusion detection. The gene fusion module uses the DRAGEN RNA spliced aligner to perform split-read analysis on supplementary (chimeric) alignments to detect potential breakpoints.	FASTQ  Optional: Custom reference Gene Annotation File (GTF, GFF, or GFF3)	<ul style="list-style-type: none"><li>• DRAGEN server</li><li>• BaseSpace Sequence Hub</li><li>• Illumina Connected Analytics</li><li>• Onboard MiSeq i100 Series, NextSeq 1000/2000 Systems, NovaSeq X Series</li></ul>
<a href="#">DRAGEN Microbial Enrichment Plus</a>	Provides secondary analysis of common Illumina infectious disease and microbiology enrichment panels (ie, VSP, RVEK, RPIP)	FASTQ, BCL	<ul style="list-style-type: none"><li>• BaseSpace Sequence Hub</li><li>• Illumina Connected Analytics</li><li>• Onboard MiSeq i100 Series</li></ul>
<a href="#">RNA-Seq visualization and statistical analysis</a>	Supports differential analysis, clustering, and data exploration plots	tsv, csv, txt, gz, FASTQ, BAM	Partek Flow
<a href="#">Illumina Connected Insights</a>	Supports streamlined interpretation and reporting from DRAGEN software for oncology research applications	VCF	Illumina Connected Insights with automated ingestion of VCF files



### Method 3: Total RNA-Seq

In addition to proteins, genes encode a vast array of nonprotein-coding elements that play an instrumental role in orchestrating how a cell is organized from a transcriptional regulation perspective. A major strength of total RNA-Seq lies in its ability to identify novel features of the transcriptome. NGS-based total RNA-Seq can sequence the whole transcriptome, including both coding and noncoding transcripts, without the limitation of probe design, delivering a high resolution, base-by-base view of coding and multiple forms of noncoding RNA activity. This provides a comprehensive picture of gene expression across the full transcriptome at a specific moment in time. With total RNA-Seq, the whole transcriptome—including both known and unknown regions—is captured.<sup>1-4</sup>

#### Relevant applications

##### Discovery of pathways associated with disease

Investigations into the transcriptomic differences between cancer samples and non-cancerous tissue have been shown to be useful in differentiating cancer subtypes, assessing the impact of mutations and identifying biomarkers and other variables. Combining total RNA-Seq with DNA sequencing and methylation analyses successfully classified meningioma tumors into distinctive molecular groups that correlated with predicted clinical outcomes, suggesting they are capable of informing medical therapies.<sup>28</sup>


##### Assessment of therapy response

Total RNA-Seq has the potential to identify genes and pathways associated with biological response or lack of response to novel drug therapies in model systems or retrospective studies of tissue samples. As part of a multiomic profiling study on tumor samples from melanoma patients undergoing anti-PD-1 immunotherapy, total RNA-Seq identified an interferon gamma (INFγ) expression pattern that, along with other biomarkers, robustly predicted successful response to immunotherapy.<sup>26</sup>



#### Step-by-step overview

Total RNA-Seq has five basic steps: RNA extraction, ribosomal RNA (rRNA) depletion (optional) and library preparation, sequencing, data analysis, and insights ([Figure 4](#)).



#### STEP 1 RNA extraction

For whole-transcriptome methods, you can use high-quality RNA or degraded samples as input. There are several options for RNA isolation kits depending on your sample type. Illumina recommends the QIAGEN RNeasy for extraction of RNA from cells and QIAGEN RNeasy FFPE for extracting RNA from FFPE tissue. Effectively evaluating RNA quality is a critical step in successful total RNA-Seq and can be achieved by measuring mean RNA fragment size before library prep.

#### STEP 2 Library preparation

Illumina Stranded Total RNA Prep with Ribo-Zero Plus enables preparation of whole-transcriptome sequencing-ready libraries from low RNA inputs and low-quality samples to support a wide range of RNA-Seq applications. The kit includes Ribo-Zero Plus for single-tube depletion of rRNA and globin RNA from multiple species to facilitate rich transcriptome analyses ([Table 9](#)).<sup>30</sup>

Table 9: Illumina Stranded Total RNA Prep with Ribo-Zero Plus specifications

Feature	Specification
RNA input amount	1 ng high-quality RNA, 10 ng FFPE RNA
Total assay time	~7 hr
Hands-on time	< 3 hr

STEP 3  
Sequencing

The sequencing system used depends on several factors, including the application, study size, throughput requirements, and more (Table 10).

Table 10: Experimental parameters for performing total RNA-Seq on different sequencing systems

System	Flow cell	Single reads per flow cell	No. of samples per flow cell <sup>a</sup>	Recommended read length
NextSeq 1000 and NextSeq 2000 Systems	P1	100M	2	2 × 100 bp
	P2	400M	8	
	P3 <sup>b</sup>	1.2B	24	
	P4 <sup>b</sup>	1.8B	36	
NovaSeq 6000 System	SP	1.6B	32	2 × 100 bp
	S1	3.2B	64	
	S2	8.2B	164	
	S4	10B	200	
NovaSeq X Series	10B	10B	200	2 × 100 bp
	25B	26B	384 <sup>c</sup>	
<div>a. Based on 50M reads per sample.</div> <div>b. P3 and P4 flow cells are only available on the NextSeq 2000 System.</div> <div>c. Based on Illumina indexes. Additional indexes can be added.</div>				

STEP 4  
Data analysis

Illumina recommends using the DRAGEN RNA pipeline for secondary analysis of whole-transcriptome libraries. The DRAGEN RNA pipeline is available on-premises with a DRAGEN Server, on the cloud-based BaseSpace Sequence Hub and Illumina Connected Analytics platforms, and onboard select sequencing systems (Table 11).

For studies that include analysis of metatranscriptomes from complex microbiome samples, such as human gastrointestinal (GI) tract, the Microbiome Metatranscriptome application on BaseSpace Sequence Hub provides taxonomic and functional profiling.

STEP 5  
Insights

For downstream tertiary analysis, output from the DRAGEN RNA pipeline can be automatically ingested into Illumina Connected Insights for variant interpretation and reporting. Results can also be transferred to Correlation Engine for gaining insights into the biological effects of gene expression changes, and to Partek Flow software for further statistical analysis and information-rich visualizations.

Table 11: Illumina total RNA-Seq analysis software

Pipeline	Application	Input	Access point
<a href="#">DRAGEN RNA</a>	<div>Offers multiple operating modes, including reference-only alignment and annotation-assisted alignment with gene fusion detection</div> <div>The gene fusion module uses the DRAGEN RNA spliced aligner to perform split-read analysis on supplementary (chimeric) alignments to detect potential breakpoints</div>	<div>FASTQ</div> <div>Optional: Custom reference Gene Annotation File (GTF, GFF, or GFF3)</div>	<ul style="list-style-type: none"><li>DRAGEN Server</li><li>BaseSpace Sequence Hub</li><li>Illumina Connected Analytics</li><li>Onboard NextSeq 1000/2000 Systems, NovaSeq X Series</li></ul>
<a href="#">Microbiome Metatranscriptomics</a>	Performs taxonomic and pathway enrichment analysis on microbiome derived RNA libraries	FASTQ	<ul style="list-style-type: none"><li>BaseSpace Sequence Hub</li></ul>
<a href="#">RNA-Seq interpretation</a>	Processes sequencing data from mRNA to estimate transcript abundance and identify differentially expressed transcripts across samples	csv, txt, xlsx files	<ul style="list-style-type: none"><li>Correlation Engine</li></ul>
<a href="#">RNA-Seq visualization and statistical analysis</a>	Supports differential analysis, clustering, and data exploration plots	tsv, csv, txt, gz, FASTQ, BAM	<ul style="list-style-type: none"><li>Partek Flow software</li></ul>
<a href="#">Illumina Connected Insights</a>	Supports streamlined interpretation and reporting from DRAGEN software for oncology research applications	VCF	<ul style="list-style-type: none"><li>Illumina Connected Insights with automated ingestion of VCF files</li></ul>



Support that never stops

Illumina strives to be the best partner possible. With a global presence, you can rely on our support to facilitate your success. Technical support is available via phone five days a week or access online support 24/7, worldwide and in multiple languages, with rapid response time near most major metropolitan areas. Illumina provides excellent product consistency, supply, and quality enabled by a mature global manufacturing infrastructure.

Trusted technology, trusted partner

As a preferred NGS platform provider, Illumina has shipped over 20,000 sequencing systems globally. Illumina NGS technology is cited in over 421,000 peer-reviewed publications—5× more than all other NGS technologies combined.<sup>31</sup> Building on decades of expertise, Illumina has a relentless commitment to innovation and building future NGS capabilities and applications.

Learn more →

[Illumina RNA sequencing](#)



Ordering information

Library preparation

Product	Catalog no.
Illumina Stranded mRNA Prep, Ligation (16 samples)	20040532
Illumina Stranded mRNA Prep, Ligation (96 samples)	20040534
Illumina RNA Prep with Enrichment, (L) Tagmentation (16 samples)	20040536
Illumina RNA Prep with Enrichment, (L) Tagmentation (96 samples)	20040537
Illumina Stranded Total RNA Prep, Ligation with Ribo-Zero Plus (16 samples)	20040525
Illumina Stranded Total RNA Prep, Ligation with Ribo-Zero Plus (96 samples)	20040529
Illumina Stranded Total RNA Prep, Ligation with Ribo-Zero Plus Microbiome (96 samples)	20072063
Illumina DNA/RNA UD Indexes Set A, Tagmentation (96 Indexes, 96 Samples) <sup>a</sup>	20091654
Illumina DNA/RNA UD Indexes Set B, Tagmentation (96 Indexes, 96 Samples) <sup>a</sup>	20091656
Illumina DNA/RNA UD Indexes Set C, Tagmentation (96 Indexes, 96 Samples) <sup>a</sup>	20091658
Illumina DNA/RNA UD Indexes Set D, Tagmentation (96 Indexes, 96 Samples) <sup>a</sup>	20091660
PhiX Control Kit, v3	FC-110-3001
a. 384 indexes total available when combining sets A, B, C, and D.	



Sequencing systems

System	Catalog no.
MiSeq i100 Plus System	20115695
NextSeq 1000 System	20038898
NextSeq 2000 System	20038897
NovaSeq 6000 System	20012850
NovaSeq X System	20084803
NovaSeq X Plus System	20084804

Data analysis

Product	Catalog no.
BaseSpace Sequence Hub Professional Annual Subscription <sup>a</sup>	20042109
BaseSpace Sequence Hub Enterprise Annual Subscription <sup>a</sup>	15066411
Illumina Connected Analytics Professional Annual Subscription	20044876
Illumina Connected Analytics Enterprise Annual Subscription	20038994
Illumina Connected Insights–Oncology Genome Equivalent Sample – VCF	20090138
Illumina Connected Insights Starter Implementation Package	20071787
Illumina Connected Insights Expanded Implementation Package	20071787 (as scoped)
Partek Flow software	Contact an Illumina sales representative
Correlation Engine	Contact an Illumina sales representative
DRAGEN Server	20051343
a. BaseSpace Sequence Hub subscriptions include complimentary iCredits for running analysis apps and data storage. Additional iCredits are available for purchase.	

Appendix

References

1. Oszolak F, Milos PM. [RNA sequencing: advances, challenges and opportunities](#). *Nat Rev Genet.* 2011;12(2):87–98. doi:10.1038/nrg2934

2. Wang Z, Gerstein M, Snyder M. [RNA-Seq: a revolutionary tool for transcriptomics](#). *Nat Rev Genet.* 2009;10:57–63. doi:10.1038/nrg2484

3. Wilhelm BT, Landry JR. [RNA-Seq—quantitative measurement of expression through massively parallel RNA-Sequencing](#). *Methods.* 2009;48:249–57. doi:10.1016/j.ymeth.2009.03.016

4. Wang Y, Mashock M, Tong Z, et al. [Changing Technologies of RNA Sequencing and Their Applications in Clinical Oncology](#). *Front Oncol.* 2020;10:447. Published 2020 Apr 9. doi:10.3389/fonc.2020.00447

5. Su Z, Labaj PP, Li S, et al. [A comprehensive assessment of RNA-Seq accuracy, reproducibility and information content by the Sequencing Quality Control Consortium](#). *Nat Biotech.* 2014;32:903–914. doi:10.1038/nbt.2957

6. Xu J, Gong B, Wu L, Thakkar S, Hong H, Tong W. [Comprehensive Assessments of RNA-seq by the SEQC Consortium: FDA-Led Efforts Advance Precision Medicine](#). *Pharmaceutics.* 2016;8(1):8. Published 2016 Mar 15. doi:10.3390/pharmaceutics8010008

7. Crist AM, Hinkle KM, Wang X, et al. [Transcriptomic analysis to identify genes associated with selective hippocampal vulnerability in Alzheimer’s disease](#). *Nat Commun.* 2021;12(1):2311. doi:10.1038/s41467-021-22399-3

8. Allen M, Wang X, Burgess JD, et al. [Conserved brain myelination networks are altered in Alzheimer’s and other neurodegenerative diseases](#). *Alzheimers Dement.* 2018;14(3):352–366. doi:10.1016/j.jalz.2017.09.012

9. Guennewig B, Lim J, Marshall L, et al. [Defining early changes in Alzheimer’s disease from RNA sequencing of brain regions differentially affected by pathology](#). *Sci Rep.* 2021;11(1):4865. doi:10.1038/s41598-021-83872-z

10. Zhang X, van Rooij JGJ, Wakabayashi Y, et al. [Genome-wide transcriptome study using deep RNA sequencing for myocardial infarction and coronary artery calcification](#). *BMC Med Genomics.* 2021;14(1):45. doi:10.1186/s12920-020- 00838-2

11. Luo X, Yin J, Dwyer D, et al. [Chamber-enriched gene expression profiles in failing human hearts with reduced ejection fraction](#). *Sci Rep.* 2021;11(1):11839. doi:10.1038/s41598-021-91214-2

12. Jain PN, Robertson M, Lasa JJ, et al. [Altered metabolic and inflammatory transcriptomics after cardiac surgery in neonates with congenital heart disease](#). *Sci Rep.* 2021;11(1):4965. doi:10.1038/s41598-021-83882-x

13. Meng Z, Chen C, Cao H, Wang J, Shen E. [Whole transcriptome sequencing reveals biologically significant RNA markers and related regulating biological pathways in cardiomyocyte hypertrophy induced by high glucose](#). *J Cell Biochem.* 2019;120(1):1018–1027. doi:10.1002/jcb.27546

14. Andersen JD, Jacobsen SB, Trudsø LC, Kampmann ML, Banner J, Morling N. [Whole genome and transcriptome sequencing of post-mortem cardiac tissues from sudden cardiac death victims identifies a gene regulatory variant in NEXN](#). *Int J Legal Med.* 2019;133(6):1699–1709. doi:10.1007/s00414-019-02127-9

15. Zhao S, Fung–Leung W-P, Bittner A, Ngo K, Liu X. [Comparison of RNA-seq and microarray in transcriptome profiling of activated T cells](#). *PLoS ONE.* 2014;9(1):e78644. doi:10.1371/journal.pone.0078644

16. Atak ZK, Gianfelici V, Hulselmans G, et al. [Comprehensive analysis of transcriptome variation uncovers known and novel driver events in T-cell acute lymphoblastic leukemia](#). *PLoS Genet.* 2013;9(12):e1003997. doi:10.1371/journal.pgen.1003997

17. Kumar-Sinha C, Kalyana-Sundaram S, Chinnaiyan AM. [Landscape of gene fusions in epithelial cancers: seq and ye shall find](#). *Genome Med.* 2015;7:129. Published 2015 Dec 18. doi:10.1186/s13073-015-0252-1

18. Docking TR, Parker JDK, Jadersten M, et al. [A clinical transcriptome approach to patient stratification and therapy selection in acute myeloid leukemia](#). *Nat Commun.* 2021;12(1):2474 doi:10.1038/s41467-021-22625-y

19. Figgett WA, Monaghan K, Ng Milica, et al. [Machine learning applied to whole-blood RNA-sequencing data uncovers distinct subsets of patients with systemic lupus erythematosus](#). *Clin Transl Immunology.* 2019;8(12):e01093

20. Lu L, Zhang H, Pang J, Hou G, Lu M, Gao X. [ERG rearrangement as a novel marker for predicting the extra-prostatic extension of clinically localised prostate cancer](#). *Oncol Lett.* 2016;11(4):2532–2538. doi:10.3892/ol.2016.4282

21. Perez-Gracia JL, Sanmamed MF, Bosch A, et al. [Strategies to design clinical studies to identify predictive biomarkers in cancer research](#). *Cancer Treat Rev.* 2017;53:79–97. doi:10.1016/j.ctrv.2016.12.005

22. Fang B, Mehran RJ, Heymach JV, Swisher SG. [Predictive biomarkers in precision medicine and drug development against lung cancer](#). *Chin J Cancer.* 2015;34(7):295–309. doi:10.1186/s40880-015-0028-4

23. Zhao X, Modur V, Carayannopoulos LN, Laterza OF. [Biomarkers in Pharmaceutical Research](#). *Clin Chem.* 2015;61(11):1343–1353. doi:10.1373/clinchem.2014.231712

For Research Use Only. Not for use in diagnostic procedures.

M-GL-00034 v2.0

22

23

M-GL-00034 v2.0

For Research Use Only. Not for use in diagnostic procedures.

24. Pon JR, Marra MA. [Driver and passenger mutations in cancer](#). *Annu Rev Pathol*. 2015;10:25-50. doi:10.1146/annurev-pathol-012414-040312
25. Illumina. Surveillance of infectious disease through wastewater sequencing Application Note. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/covidseq-wastewater-epidemiology-app-note-m-gl-00429.pdf](#). Published 2022. Accessed April 24, 2024.
26. Illumina. Rapid detection of respiratory pathogens using the MiniSeq System Application Note. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/respiratory-pathogen-panel-miniseq-app-note-m-gl-01508.pdf](#). Published 2023. Accessed April 24, 2024.
27. Illumina. Illumina RNA Prep with Enrichment Data Sheet. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-rna-prep-enrichment-data-sheet-m-gl-02145.pdf](#). Published 2023. Accessed April 24, 2024.
28. Nassiri F, Liu J, Patil V, et al. [A clinically applicable integrative molecular classification of meningiomas](#). *Nature*. 2021;597(7874):119-125. doi:10.1038/s41586-021-03850-3
29. Newell F, Pires da Silva I, Johansson PA, et al. [Multiomic profiling of checkpoint inhibitor-treated melanoma: Identifying predictors of response and resistance, and markers of biological discordance](#). *Cancer Cell*. 2022;40(1):88-102.e7. doi:10.1016/j.ccell.2021.11.012.
30. Illumina. Illumina Stranded Total RNA Prep, Ligation with Ribo-Zero Plus Data Sheet. [illumina.com/content/dam/illumina/gcs/assembledassets/marketing-literature/illumina-stranded-total-rna-prep-data-sheet-m-gl-02148.pdf](#). Published 2020. Accessed August 20, 2024.
31. Data calculations on file. Illumina, Inc. 2022.



1.800.809.4566 toll-free (US) | +1.858.202.4566 tel  
techsupport@illumina.com | www.illumina.com

© 2025 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see [www.illumina.com/company/legal.html](#).  
M-GL-00034 v2.0