

MiniSeq[™] Sequencing System

The power of proven Illumina sequencing in an approachable and easy-to-use research tool.

Highlights

- Accessible Illumina sequencing
 - Affordable to acquire and cost-effective to run, even with low numbers of samples
- Push-button operation and simple data analysis
 - Walk away library-to-results solution with onboard data
- Highly flexible to fit research demands Supports a broad range of DNA and RNA sequencing applications for examining single genes to entire pathways
- End-to-end support Illumina scientists and engineers provide installation, training, and support, from assay design through data analysis



Introduction

The MiniSeq System (Figure 1) delivers the quality and reliability of Illumina next-generation sequencing (NGS) technology in a powerful, accessible benchtop sequencer with a small footprint. It enables researchers to take control of their sequencing projects. This small, robust system turns a broad range of NGS methods into approachable, easy-to-use research tools. With the MiniSeq System, there is no need to wait to batch samples for sequencing on a high-throughput instrument; researchers can sequence on demand. It circumvents the iterative, time-consuming testing of Sanger sequencing and gPCR to allow for interrogation of individual genes to entire pathways with full-gene coverage. Laboratories of any size can perform a range of sequencing methods to deliver results and advance their research.

Figure 1: The MiniSeq System — By harnessing advances in SBS chemistry and simple, streamlined workflows, the MiniSeq System delivers a library-to-results solution that is powerful and easy to use.

Powerful sequencing made simple

The MiniSeq System features a simple, integrated, library-toresults workflow that enables sequencing of both DNA and RNA with minimal hands-on time (Figure 2). It is ideal for targeted research applications such as cancer sequencing and gene expression profiling. Onboard, touch-screen data analysis with a simple, intuitive user interface eliminates the need for specialized equipment or bioinformatics expertise. Illumina scientists are available at every point along the way with support and guidance, enabling researchers to focus on making the next breakthrough discovery.




Figure 2: MiniSeq System sequencing workflow—The MiniSeq System offers a simple, integrated workflow from library preparation to onboard data analysis. Workflow times will vary by experiment and assay type. Details shown are for a sequencing run using the AmpliSeg™ for Illumina Sequencing Solution and a read length of 2 × 150 bp.

Table 1: Flexibility for multiple applications

	High-output reagent kit		Mid-output reagent kit	
Application	No. of samples	Run time ^a	No. of samples	Runtime
Targeted DNA amplicon sequencing				
207 amplicons	96	24 hours	32	17 hours
500× coverage	96	24 Hours		17 Hours
$2 \times 150 \text{ bp}$				
Targeted expression profiling				
65 targets	384	7 hours	123	6 hours
1 × 50 bp				
Enrichment panel				
1 Mb region	00	13 hours	7	10
100× coverage	23			12 hours
2×75 bp				
Viral Pathogen Panel				
1M reads/sample	20	< 5 hours	N/A ^b	N/A ^b
1 × 100 bp (Rapid Kit)				
microRNA sequencing				
5M reads/sample	5	4 hours	2	4 hours
1 × 36 bp				
Small whole-genome sequencing				
5 Mb genome	50	O4 hauwa	16	17 ho:::::
30× coverage	50	24 hours	16	17 hours
2 × 150 bp				

a. Run times are without indexes.

Streamlined sequencing workflow

The MiniSeg System provides an intuitive user interface and loadand-go operation, making it easy to learn and easy to use. It integrates clonal amplification, sequencing, and data analysis into a single instrument, eliminating the need to purchase and operate specialized, ancillary equipment. After library preparation using a simple, streamlined Illumina library prep kit, libraries are loaded into the MiniSeg System where sequencing is automated. It takes less than five minutes to load and set up a run on the MiniSeq System. Runs are complete in less than a day, and data analysis is performed onboard the instrument or in BaseSpace[™] Sequence Hub - the Illumina genomic computing environment. A suite of data analysis tools and a growing list of third-party BaseSpace Applications (Apps) empower researchers to perform their own informatics analysis easily.

By employing industry-leading Illumina sequencing by synthesis (SBS) chemistry and file format conventions, the MiniSeq System offers customers access to a broad ecosystem of established protocols, workflows, data sets, and data analysis tools.

Supports a wide range of applications

The MiniSeq System combines industry-leading Illumina NGS technology with a broad range of library preparation and data analysis solutions to deliver robust NGS tools in a simple, intuitive user experience. It offers cross-method flexibility, enabling easy transition between sequencing projects for both DNA and RNA applications. Demonstrated and optimized workflows are available for small RNA discovery, targeted resequencing, targeted RNA sequencing, and profiling of solid and hematological tumors (Table 1).

The MiniSeq System delivers a < 1-day turnaround for numerous sequencing methods. The output of the system allows researchers to sequence a broad range of samples per run:

- 1-96 targeted panel samples
- 1–384 gene expression profiling samples
- 1-12 small RNA (miRNA) profiling samples
- 1-20 viral pathogen RNA enrichment samples

The MiniSeq System is supported by the full suite of Illumina library preparation solutions, enabling library compatibility across the Illumina sequencing portfolio. This allows researchers to scale up studies easily to the higher throughput NextSeq™ Series of Sequencing Systems or perform follow-up studies on the MiSeq[™] Series of Sequencing Systems (Figure 3).

Industry-leading SBS chemistry delivers high accuracy

At the core of the MiniSeq System is industry-leading Illumina SBS chemistry, the most widely adopted NGS technology worldwide.¹ This proprietary reversible terminator-based method enables the massively parallel sequencing of millions of DNA fragments, detecting single bases as they are incorporated into growing DNA strands. The method significantly reduces errors and missed calls associated with strings of repeated nucleotides (homopolymers). The low cost-per-base allows deeper sequencing for more sensitivity and greater accuracy (Table 2).

b. N/A = not applicable

Table 2: MiniSeq System performance parameters

Flow cell configuration ^a	Read length (cycles)	Output (Gb)	Run time ^b	Data quality ^c
High-output kit	300	~ 7.5	~ 24 hours	Q30 > 80%
Up to 25M single reads	150	~ 4	~ 13 hours	Q30 > 85%
Up to 50M paired-end reads	75	~ 2	~ 7 hours	Q30 > 85%
Rapid kit	100	~ 2	< 5 hours	Q30 > 85%
Up to 20M single reads	100			
Mid-output kit				
Up to 8M single reads	300	~ 2.5	~ 17 hours	Q30 > 80%
Up to 16M paired-end reads				

- a. Actual performance parameters may vary based on sample type sample quality, and clusters passing filter.
- b. Times include cluster generation, sequencing, and base calling with quality scores on a MiniSeq System.
- c. The percentage of bases > Q30 is averaged over the entire run.

Push-button data analysis and streamlined bioinformatics

The MiniSeq System features onboard data analysis in an intuitive user interface. The instrument computer processes base calls and quality scores generated during the sequencing run. Researchers have several options for data analysis.

Local Run Manager software is a multifunctional, integrated onboard solution. Local Run Manager not only allows users to create a sequencing run, monitor status, and view results, but also analyze data. It is easily accessed through a web browser and integrates with the instrument control software. Samples to be sequenced and analysis input files are recorded, and onboard data analysis is automatically performed upon completion of the sequencing run. This produces alignment information, structural variants, expression analysis, small RNA analysis, and more for each sample based on the user-specified analysis workflow.

Also, sequencing data can be run through a wide range of opensource or commercial pipelines developed for Illumina data, or instantly transferred, analyzed, archived, and shared securely with BaseSpace Mequence Hub. BaseSpace Sequence Hub is the only cloud ecosystem that offers direct instrument integration, enabling automatic encrypted data flow directly from the instrument into the cloud ecosystem for analysis, storage, sharing, and other forms of data management. Additionally, BaseSpace Sequence Hub users can monitor the status of their runs through the cloud portal or through the iOS app for BaseSpace.



Figure 3: Illumina NGS portfolio of benchtop sequencers—Illumina NGS systems offer solutions for virtually every application, sample type, and sequencing scale. Each delivers high data quality and accuracy with flexible throughput and simple, streamlined workflows. Data can be easily compared, exchanged, and analyzed in BaseSpace Sequence Hub.

Summary

The MiniSeq System is a small, robust benchtop sequencer that enables NGS to become an everyday tool in laboratories worldwide. Incorporating advances in SBS chemistry, the flexible MiniSeq System features push-button operation and streamlined library-to-results workflows that allow researchers to perform popular NGS applications. Its price point and cost-effective operation, even for low numbers of samples, makes the power of proven Illumina sequencing more accessible than ever.

Ordering Information

Order the MiniSeg System and reagents at www.illumina.com

System	Catalog no.	
MiniSeq Sequencing System	SY-420-1001	
Sequencing reagent kits	Catalog no.	
MiniSeq High Output Kit (75 cycles)	FC-420-1001	
MiniSeq High Output Kit (150 cycles)	FC-420-1002	
MiniSeq High Output Kit (300 cycles)	FC-420-1003	
MiniSeq Rapid Kit (100 cycles)	20044338	
MiniSeq Mid Output Kit (300 cycles)	FC-420-1004	

MiniSeq System specifications

Parameter	Specifications	
Instrument configuration	RFID tracking for consumables	
Instrument control computer (internal) ^a	Base unit: Intel Core i7-4700EQ 2.4 GHz CPU Memory: 16 GB DDR3L RAM Hard drive: 1 Tb Operating system: Windows 10 embedded standard	
Operating environment	Temperature: 19°C to 25°C (22°C ± 3°C) Humidity: noncondensing 20%–80% relative humidity Altitude: less than 2000 m (6500 ft) Air quality: pollution degree rating of II, air particulate cleanliness levels per ISO9 (ordinary room air) or better Ventilation: up to 2048 BTU/hr @ 600 W For indoor use only	
Light emitting diode (LED)	515 nm, 650 nm	
Dimensions	W×D×H: 45.6 cm × 48 cm × 51.8 cm (18.0 in × 18.1 in × 20.4 in) Weight: 45 kg (99 lbs) Crated weight: 56.5 kg (125 lbs)	
Power requirements	100-120 volts AC — A 15 Amp grounded 220-240 volts AC — A 10 Amp grounded	
Radio frequency identifier (RFID)	Frequency: 13.56 MHz Power: supply 3.3 volts DC ± 5%, current 120 mA, RF output power 200 mW	
Product safety and compliance	NRTL certified IEC 61010-1 CE marked to the Low Voltage Directive 2006/95/EC FCC/IC approved	

Learn more

To learn more about the MiniSeq System, visit www.illumina.com/systems/sequencing-platforms/miniseq.html

To learn more about sequencing applications on the MiniSeq System, visit MiniSeq System Literature to access relevant application notes

References

1. Data calculations on file. Illumina, Inc., 2017.

Maximize Performance and Productivity with Illumina Services, Training, and Consulting

Whether immediate help is needed during an instrument run, or in-depth consultations are required for sophisticated workflows, Illumina can help. Illumina service and support teams provide a full suite of expedient, customized solutions from initial trainings, to instrument support, and ongoing NGS consultations. Illumina Professional Care Services Packs

Illumina offers Professional Care Services Packs, allotments of points that can be redeemed for discounted Illumina Professional Services. Benefits of Services Packs include:

- One-time investment no need for additional, postsale expenditures
- Risk mitigation bank points for unanticipated future services
- Savings cost-effective versus a la carte pricing

Professional Services

Product Care Services

- Tiered Instrument Service Plans + Add-On Services
- Instrument Compliance Services
- Instrument On-Demand Services

Illumina University Training

- Instructor-Led Training at Your Chosen Facility
- Instructor-Led Training at an Illumina Training Center
- Online Courses
- Webinars

Illumina Consulting

- Proof-of-Concept Services for instrument and library preparation testing
- Concierge Services for design assistance and product optimization
- IT and Bioinformatics Hourly Consulting for personalized assistance

For more on Illumina support offerings, visit:

www.illumina.com/services/instrument-services-training.html

