


Figure 1: Detection of Changes in DNA Copy Number Using Human-1 and Humanhap300 Genotyping Beadchips



Single copy deletions, monoallelic duplications, and amplifications are modeled on the Human-1 BeadChip (109K; A-D) and the HumanHap300 BeadChip (317K; E) using cell lines with one to four X chromosomes. All plots are shown juxtaposed with normal genomic profiles from chromosome 10 and were created with Illumina BeadStudio software.

(A) For XY, the presence of a single X chromosome is shown as a decrease in the log R ratio. In the Allele Frequency (AF) plot, the heterozygous state completely collapses to the homozygous axis.

(B) For XX, the presence of the expected two copies of the X chromosome shows no deflection in the log R ratio. In the heterozygotes are clustered around +0.5.

(C) For XXX, the log R ratio increases. The heterozygous state splits into two clusters representing a 2:1 and 1:2 ratio.

(D) For XXXX, the log R ratio increases even further. The heterozygous state is divided into three populations with allelic ratios of 3:1, 2:2, and 1:3, respectively.

(E) The response of the log R ratio for both the X chromosome and chromosome 10 for each X-copy number cell line on both the Human-1 and HumanHap300 Genotyping BeadChips. The log R ratio increases with increasing X-copy number of the X chromosome but not for chromosome 10. The corresponding standard deviation is shown. Note the similarities between Infinium I and Infinium II Assay-derived results. This data was generated with a 10-SNP moving average.

Examples of Allelic Imbalance in Cell Lines


The ability of the HumanHap550 Genotyping BeadChip to detect allelic imbalances in the cancer cell lines HL60 (promyelocytic leukemia) and MCF7 (breast adenocarcinoma) is shown in Figure 2. All of the plots shown were generated with the Illumina Genome Viewer in BeadStudio. In HL60, deletions of varying size appear on Chromosome 9: a ~20Mb deletion on the p arm, and a ~1.8Mb deletion on the q arm (Figure 2A). These hemizygous deletions (from two copies to one copy) are manifest as a downward deflection in the log R ratio and a loss of heterozygotes in the AF. Other examples of aberrations in HL60 include monoallelic duplication of Chromosome 18 as indicated by an increase in the log R ratio and the split of the heterozygotes into two states: one at 0.67 (2:1 ratio) AF and another at 0.33 (1:2 ratio) AF (Figure 2B). There are also several small amplifications of approximately 260kb in HL60 on Chromosome 8 with 1:9 and 9:1 allelic ratios (Figure 2C). All of these aberrations have previously been identified using SKY karyotyping⁸. Adjacent homozygous and hemizygous deletions are revealed in MCF7 (Figure 2D). The gene annotation

provided in the Illumina BeadStudio software allows users to quickly analyze affected gene regions. In this case, both copies of *MTAP* are homozygously deleted and *CDKN2A* (p16) and *CDKN2B* (p15) are hemizygously deleted.

Summary

Illumina's Infinium Assays and Whole-Genome Genotyping BeadChips offer genome-wide coverage at an unprecedented resolution. The combined measurement of allelic ratios and normalized intensities provides enhanced detection of aberrations while facilitating identification of copy-neutral genetic anomalies such as uniparental disomy (UPD) and mitotic recombination^{9,10}. Whole-Genome Genotyping can also yield allelic information on deletions, duplications, and amplifications, which have implications in cancer therapeutics. The combination of proven assays, high-density arrays, and integrated software enables the analysis of LOH and copy number changes in both single and paired (i.e., matched) samples with high precision.

Figure 2: Detection of Allelic Imbalances in the Cancer Cell Lines HL60 and MCF7



Genomic profiles of cancer cell lines on the HumanHap550 Genotyping BeadChip generated with Illumina BeadStudio software.

(A) Two regions in HL60 that exhibit LOH are shown by a downward deflection in the log R ratio and loss of heterozygotes in the AF profile.

(B) The entire length of Chromosome 18 in HL60 exhibits a monoallelic duplication (trisomy), as evidenced by the elevated log R ratio and split heterozygous cluster.

(C) Several amplifications in HL60 are present on Chromosome 8.

(D) High-resolution analysis of MCF7 showing a 150kb homozygous and 83kb hemizygous deletion in a region on Chromosome 9 containing *MTAP*, *CDKN2A* (p16), and *CDKN2B* (p15). For all log R ratio profiles, the blue line indicates a 100kb moving median for the HumanHap550 Genotyping BeadChip.

