

# Sequencing Analysis Viewer Software Release Notes

## SAV 1.9.1

*For all sequencing systems*

**March 5, 2015**

FOR RESEARCH USE ONLY

© 2015 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CPro, CytoChip, DesignStudio, Epicentre, GAIIX, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, ForenSeq, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, Nextera, NextBio, NextSeq, Powered by Illumina, SeqMonitor, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

## Introduction

These Release Notes detail the key changes to software components for the Sequencing Analysis Viewer (SAV) Software since the package containing v1.8.46.

This SAV version requires Windows .NET Framework 4.5.1, systems using Windows .NET Framework 4 should remain on SAV 1.8.37 or prior. See this link for more information.

- Windows .NET Framework 4.5.1 Offline installer: <http://www.microsoft.com/en-us/download/details.aspx?id=40779>

If you are upgrading from a version before SAV v1.8.46, review the release notes for SAV v1.8.46 for a list of features and bug fixes introduced in that version.

The software package includes:

- Sequencing Analysis Viewer (SAV) Software v1.9.1

## I. Sequencing Analysis Viewer 1.9.1

### NEW FEATURES:

- None

### IMPROVEMENTS:

- None

### DEFECT REPAIRS:

- Fixed exception handling of issue if user tries to access the clipboard at the same time to prevent SAV crash.
- Fixed table and graph selection/refresh bugs.

### KNOWN ISSUES:

- Windows Control Panel, Programs and Features can show entries for previously installed version after updating.