

Decode File Client Release Notes

January 2024

Introduction

Decode File Client (DFC) provides a way to search, download, and view DMAP files for use with Illumina® microarray scanners. DMAP files are the unique map files that provide the location of probes for each Infinium™ BeadChip and are required for scanning. The software requires internet access and permission to write files to a specific folder upon download.

DMAP files are available 24 to 48 hours after shipment of all Infinium BeadChip orders. The files remain accessible until 2 years after the date of manufacture.

These Release Notes detail the key features and changes to software components for the release of Decode File Client. For information on how to use the system, see the [Decode File Client Software Guide](#).

v3.0.4 NEW FEATURES IN DETAIL

- Updated installer executable to improve installation experience in certain Windows environments.
 - Successful installations of v3.0.3 do not require v3.0.4 update.

v3.0.3 NEW FEATURES IN DETAIL

- Transport Layer Security (TLS) upgraded from version 1.1 to version 1.2.
 - Communication using TLS v1.1 will no longer be supported by Amazon Web Services (AWS).
- Installer executable upgraded to detect previous version 3.0.2 and provide upgrade so deleting prior versions of the software is not required after upgrading.

KNOWN ISSUES AND LIMITATIONS

- Maximum of 300 BeadChips can be downloaded without error.