

# Symbol Key

## FOR IN VITRO DIAGNOSTIC USE

This is a general legend to assist comprehension of symbols that may appear in the labeling of an Illumina *in vitro* diagnostic (IVD) product.

	Contains sufficient for <n> tests
	<i>In vitro</i> diagnostic medical device
	Authorized Representative in the European Community
	European Conformity
	Part Number
	Serial Number
	Temperature limit
	Batch code
	Manufacturer
	Date of manufacture
	Use-by date
	Consult instructions for use
	Used for medical devices that require protection from light or heat
	Reagents are to be used only in the Post-Amplification area of the laboratory
	Reagents are to be used only in the Pre-Amplification area of the laboratory
	Caution – Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner



Health Hazard Symbol - For an explanation of the hazard symbols and the accompanying hazard statement please refer to component Safety Data Sheet (SDS) which can be downloaded from <http://support.illumina.com/sds.html>



Hazard Warning Symbol - For an explanation of the hazard symbols and the accompanying hazard statement please refer to component Safety Data Sheet (SDS) which can be downloaded from <http://support.illumina.com/sds.html>


## Copyright and Trademarks

© 2019 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see [www.illumina.com/company/legal.html](http://www.illumina.com/company/legal.html).

## Contact Information



Illumina  
 5200 Illumina Way  
 San Diego, California  
 92122 U.S.A.  
 +1.800.809.ILMN (4566)  
 +1.858.202.4566 (outside  
 North America)  
 techsupport@illumina.com  
 www.illumina.com



Illumina Cambridge Limited  
 Chesterford Research Park,  
 Little Chesterford  
 Saffron Walden, CB10 1XL  
 UNITED KINGDOM

**Australian Sponsor**  
 Illumina Australia Pty Ltd  
 1 International Court  
 Scoresby, Victoria, 3179  
 Australia