

Correlation Engine

Customer Release Notes

Version 2.2.0

August 28, 2023

Introduction

These Release Notes detail the latest release of Correlation Engine, including known issues.

Correlation Engine provides an interface with tools and content for exploratory, preclinical, and early translational research in drug discovery and biotechnology research. The software is a multidimensional integration platform designed to search and correlate global collections of heterogeneous, large-scale experimental data across diverse platforms, organisms, and therapeutic areas.

NEW FEATURES:

- Swagger integration for Application Programming Interfaces (APIs)
- Updated FAQ
- Updated Terms of Service
- Pathway Viewer and Pathway Studio links deprecated
- mSigDB collection 8 cell types addition
 - Note: New imports will show results with the new mSigDB collection. Older imports will show incomplete results while re-scoring against the new biogroups takes place over the next several weeks.

FIXED ISSUES:

- Fixed an issue where "Edit description" from Action drop down in Curated Studies was disabled
- Fixed an issue where data filter did not return cell line results when filtering for copy number or somatic mutations in Cell Lines – Body Atlas

KNOWN ISSUES:

- No reported issues

DATA CONTENT UPDATE:

	Added since last release	Total public content
Studies	515	25,623
Biosets	13,071	186,929

Release History

Revision	Release Reference	Originator	Description of Change
00	1090760	Joe Delaney	Initial Release