

# **iScan Control Software (ICS)**

## **Customer Release Notes**

**iCS Version 4.3.0**

*For iScan systems*

**February 18, 2025**

## Introduction

These Release Notes detail the key changes to software components for the iScan Control Software since the package containing ICS 4.2.1

This is an optional update meant for customers who run EX format beadchips.

If you are upgrading from a version prior to ICS 4.3.0 please review the release notes on Illumina's website for a list of features and bug fixes introduced in those versions.

ICS 4.3.0 introduces a post-processing image sharpening method to improve consistency of data quality when scanning EX format BeadChips

The software package includes:

- iScan Control Software v4.3
- Universal Copy Service 2.7.3

## I. iScan Control Software v4.3

### NEW FEATURES:

- Implements post-processing image sharpening method for EX format BeadChips
- If iScan is configured to send data to ICA, ICA v2 will be used

### FIXED ISSUES:

- N/A

### KNOWN ISSUES:

- N/A

## Release History

[This is a **mandatory** section. It should be the very last section of a document.]

Revision	Release Reference	Originator	Description of Change
00	CN 1119024	Rey Sequerra	Initial Release