

cBot Software Release Notes

cBot v3.0.56

***For cBot/cBot2 instruments used with HiSeq 2500,
HiSeq 3000/4000, and HiSeq X systems***

June 16, 2017

FOR RESEARCH USE ONLY

Page **1** of **2**

© 2017 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CPro, CytoChip, DesignStudio, Epicentre, GAIIX, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, ForenSeq, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, Nextera, NextBio, NextSeq, Powered by Illumina, SeqMonitor, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

Introduction

These release notes describe key changes to software components for cBot v3.0.56 since the release of cBot v3.0.46.

The software package includes:

- cBot v3.0.56 installer

I. cBot 3.0.56

NEW FEATURES:

- Added SR and PE custom primer recipes for HiSeq3000/4000
 - *HiSeq_3000_4000_PE_Exclusion_Amp_v1.0_custom_primer.xml*
 - *HiSeq_3000_4000_SR_Exclusion_Amp_v1.0_custom_primer.xml*
- Enabled Illumina Advantage for all reagent kits (TG)

KNOWN ISSUES:

- None