

Software Release Notes

Local Run Manager

RNA Fusion Module

TBD

Introduction

These Release Notes detail notable items for the Local Run Manager RNA Fusion Module release.

The Local Run Manager software is an integrated solution for recording samples for a run, specifying run parameters, monitoring run status, performing data analysis, The Local Run Manager RNA Fusion analysis module aligns reads against a specified reference using the STAR aligner, and then detects gene fusions using Manta. This workflow is designed specifically for RNA libraries prepared with the TruSight RNA.

I. Local Run Manager RNA Fusion Module version 1.0.0.351

NEW FEATURES:

- Initial Release

DEFECT REPAIRS:

- None

KNOWN ISSUES:

- None